



# INTRUSAFE

CARE + MAINTENANCE PRODUCTS

---

## Instrument Identification Tape

The Instrument Identification Tapes are intended to identify surgical instruments through processing and use. The Tape material is vinyl with acrylic adhesive.

They have been evaluated for the following sterilization methods:

- Conventional Steam Autoclave Cycle for 121° C for 30 minutes
- Steam Autoclave Cycle at 132° for 3 minutes
- Conventional Hospital Ethylene Oxide Cycle

The Tape is gas and steam permeable.