

InstruSafe® Trays with Wrap + Containers



FDA 510(K) Clearance Quick Reference Chart

InstruSafe Instrument Protection Trays have been cleared for each sterilization cycle listed below when using wrap as a sterile barrier and / or as indicated, when using with rigid containers.

STERILIZATION CYCLE	WRAP	AESCULAP® RIGID CONTAINER	GENESIS® RIGID CONTAINER
Steam 4 minute	X	X	X
Ethylene Oxide	X		X
STERRAD® 100S Standard	X	X	
STERRAD® 100NX® Standard	X	X	
STERRAD® 100NX® Express	X		
STERRAD® 100NX® Flex	X	X	
STERRAD® NX®	X	X	
STERIS® AMSCO® V-PRO®, VPRO® 1 + V-PRO® maX	X	X	

For compatible rigid container size information, see Appendix A in the InstruSafe IFU. All information is subject to change, please reference www.instrusafe.com for the most up to date information.

Genesis® is a registered trademark of CareFusion Corporation. Aesculap® is a registered trademark of B. Braun Company. STERRAD®, 100S®, 100NX® are registered trademarks of Advanced Sterilization Products (ASP) Division of Ethicon, Inc. STERIS®, V-PRO®, V-PRO® 1 and V-PRO® maX are registered trademarks of STERIS Corporation.

Summit Medical is an FDA registered facility, ISO 13485 : 2012 certified and CMDCAS certified.
©2018 Summit Medical. InstruSafe® is a registered trademark of Summit Medical. All rights reserved.

QUAL-4088 | RevD | AUG2018