



INTRUSAFE

INSTRUMENT PROTECTION TRAYS

Your Trusted Leader for Instrument Longevity.



About Innovia Medical

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com

THE INSTRUSAFE ADVANTAGE	
360° OF PROTECTION	4
INDUSTRY DISADVANTAGES	6
FEATURES + BENEFITS	8
ORDER PROCESS: THE INSTRUSAFE WAY	10
STERILIZATION BARRIERS + FOOT OPTIONS	12
TRAYS BY PROCEDURE	
ARTHROSCOPY	14
BARIATRIC	16
CARDIOVASCULAR	18
DENTAL	20
ENT	22
GENERAL	34
GYNECOLOGY	38
NEUROLOGY/SPINE	40
OPHTHALMOLOGY	46
ORTHOPEDIC	50
PLASTICS	52
RETRACTORS	54
ROBOTICS	58
SCOPES	70

TRANSPORT CONTAINERS	
XL TRANSPORT CONTAINERS	80
HEAVY DUTY TRANSPORT CONTAINERS	82
TRANSPORT CONTAINERS	83
ACCESSORIES	
SILICONE SMALL PARTS BOXES	84
PLAQUE OPTIONS	85
FINGER MATS	86
GUARDIAN TRAY	86
INSTRU-ZYME	87
ROUNDERS POINT GUARDS	87

RESOURCES

POINT CAMERA TO SEE VALIDATIONS

POINT CAMERA TO READ IFU

SUBMIT PURCHASE ORDERS TO:

US + INTERNATIONAL ORDERS
Summit Medical
 815 Vikings Parkway, Ste 100
 St. Paul, MN 55121
 USA

1-888-229-2875 | +1 651-789-3939
www.innoviamedical.com/instrusafe

UK ORDERS
DTR Medical
 17 Clarion Court, Enterprise Park
 Swansea SA6 8RF
 UK

+44 (0) 1792 797910
www.dtrmedical.com

The InstruSafe® Advantage

360° of Protection

The use of complex and delicate instrument sets is a new reality in which sterile processing and the operating room live. Intricacy of surgical instrumentation continues to increase, and so does the cleaning and sterilization process. Now, healthcare facilities face the challenge of reprocessing instruments safely and efficiently while also meeting the demand for fast turnaround. Keeping these sets organized and protected upon arrival in the operating room creates a challenge for all involved. Relieve yourself of these daily obstacles with InstruSafe Trays—an everyday solution for the intricacies of your instrument sets.

We have spent over 30 years understanding the use cycle and improving outcomes. With the use cycle in mind, we have created a product that eliminates common frustrations associated with instrument sets. The products in this catalog are a compilation of our hard work. As you journey through the pages, let us help you solve your everyday instrument problems.

Summit Medical's most recent FDA 510(k) clearances and product advancements are a testament of our dedication to our customers' satisfaction. We wouldn't have continuing success without the healthcare community we serve. Focus, determination and your support will continue our impact on patient care for decades to come.



Designed to organize and protect from the point of operating room use to decontamination, our solutions benefit everyone who handles our products.

FDA 510(k) Clearance

Designed to ensure optimal sterilization, InstruSafe Trays have undergone rigorous testing before receiving multiple FDA 510(k) clearances. Cleared for various sterilization cycles, these reliable trays give hospitals and healthcare facilities even more options when it comes to organizing, sterilizing, and protecting their instruments.

For the most up-to-date list of cleared cycles, visit www.instrusafe.com/validations.



POINT CAMERA TO SEE VALIDATIONS



DECONTAMINATION

Instruments are individually protected going into the wash cycle and the trays are designed with large perforations, allowing maximum detergent and water flow.

TRAY PREPARATION

Each tray is designed specifically for the instruments within it, providing individual locations for each one. This protects them from harming the technician, from damaging themselves and from puncturing or jeopardizing the integrity of the sterile barrier.

STERILIZATION

The large amount of perforations in our trays provide superior sterilant penetration. Many of our silicone parts are also designed with minimal contact points and additional ventilation in each slot that match our tray holes to improve sterilant penetration to each individual instrument.

TRANSPORTATION

Soft and durable silicone holders in the base and cover secure instruments to prevent shifting and provide 360 degrees of protection. InstruSafe Trays help to avoid unnecessary damage and to present an organized instrument set every time.

OPERATING ROOM USE

Each instrument is easily identifiable and placed in a singular slot for surgical presentation and tracking at the end of a procedure.

Look Familiar?

Unfortunately, this is what most sterile processing departments have come to know. The familiarity and acceptance of these messy, hazardous baskets makes us cringe! Why, you ask? Because the minimal organization and protection of surgical instruments practiced in many facilities presents numerous risks to the individuals that come in contact with them, as well as to the instruments themselves.

DISADVANTAGES OF MESSY BASKETS

POTENTIAL TO COMPROMISE STERILE BARRIER BY POKING THROUGH WRAP

DANGEROUS TO TRANSPORT FROM THE O.R. TO THE SPD

MINIMAL INSTRUMENT PROTECTION AND ORGANIZATION

DIFFICULT TO COUNT AND TRACK INSTRUMENTS

BASKETS NEED FREQUENT REPLACEMENT

ADDS TIME TO SURGICAL SET-UP

BROKEN INSTRUMENTS IF DROPPED

HIGH RISK FOR FINGER POKES

POOR SURGICAL PRESENTATION

HARD TO FIND WHAT YOU NEED

MISSING/LOST INSTRUMENTS

UNNECESSARY REPROCESSING

DAMAGED INSTRUMENTS

WASTE OF SPACE




EXPECT MORE.


For these reasons and many more, InstruSafe Trays offer custom solutions to meet each customer's specific needs. We think it's time for everyone to expect more from their trays.

Absolute Protection

Trusted Solutions + Superior Results




QUARTER TURN LATCH
Secures the tray cover to the tray base, providing visual indication of the open/close position. Our off-set latch hole aligns the cover in the correct position ensuring the silicone holders match up with the instruments below.



RETRACTABLE HANDLES
Our multi-functional space-saving design fits into rigid containers making it effortless to remove the tray.



TF FOOT
Our Tall Foot fits into rigid containers providing clearance for the filter in the base.




R FOOT
Rounders Corner Guards prevent tears in sterilization wrap during storage and transportation.


FOOT OPTIONS
Sterile barrier and storage needs determine the foot type built into a tray.






PERFORATIONS (GRID PATTERN)
Each tray contains a large quantity of perforations that provide maximum sterilant penetration and ventilation.



VENTILATION
Our unique silicone design has channels that align with our perforations allowing detergent and sterilant to flow through the holder and access the entire instrument.



SILICONE HOLDERS
InstruSafe silicone holders are soft yet durable, with a history of longevity. Our variety of holders accommodate instruments large or small, allowing each instrument to have its own location within the tray. Holders in the tray base and cover provide 360 degrees of protection reducing instrument repair and replacement while delivering instruments to the O.R. in their optimal working condition.



REDUCED CONTACT POINTS
Each silicone holder is individually designed with ridges and channels built to reduce contact with the instrument while securely holding the instrument in place.

5 Simple Steps

If one of our preconfigured tray designs does not fit your needs let our product experts design a custom tray with our effortless process below.




1 ASSESS

An InstruSafe Representative conducts a free, on-site assessment of your instrument set.



2 GATHER

Your instrument set specifics and inputs are gathered by our team of experts.



3 DESIGN

Using your inputs, our product specialists build a tray design with you and your instrument set in mind.



4 EVALUATE

You be the judge! Evaluate the new tray design while it's loaded with your instrument set. Let us know your thoughts and we will make any and all adjustments necessary!



5 COMPLETE

It's all yours! Once complete, each custom tray is assigned its own unique part number, making it simple and easy for you to reorder in future.

HOW TO USE THE CATALOG


The colored headers identify which procedure the instrument set you're looking at is used for. The trays on each page are differentiated by preconfigured or custom configured with the symbols shown in the key below. Some trays have several layers including a cover, shelves, baskets, and a base. The custom configured tray examples included in the catalog are a small representation of the trays our local representatives have built for facilities around the world.

CUSTOM + PRECONFIGURED SYMBOLS




THIS SYMBOL IDENTIFIES A PRECONFIGURED TRAY

Order by part number from your local InstruSafe representative or call us today.



THIS SYMBOL IDENTIFIES A CUSTOM CONFIGURED TRAY EXAMPLE

Call us today to schedule your free assessment.



Sterilization Barriers

InstruSafe Instrument Protection Trays are FDA cleared for several sterilization cycles using wrap and/or rigid containers.

For the most up-to-date list of FDA 510(k) cleared cycles, visit www.instrusafe.com.

Foot Options

R Foot: Rounders Corner Guards—to be used with wrap


TF Foot: Tall Foot fits into rigid containers

Rounders® Corner Guards

Rounded, medical grade corner guards affix to InstruSafe Tray corners to reduce potential damage to the sterile wrap barrier. Rounders are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to all InstruSafe Trays.

Tall Foot

Designed to affix to trays that are intended to be sterilized within a rigid container. Tall Feet are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to most InstruSafe Trays.



ROUNDERS CORNER GUARDS

Small Rounders	
Kit of Four	IN-7334
Large Rounders	
Kit of Four	IN-7337

TF FOOT KIT

Kit of Four	IN-7338
-------------	---------

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit all InstruSafe Trays
- Rounded, medical grade silicone corner guard
- Help avoid compromised sterility in the O.R.
- Environmentally and economically efficient
- Impact resistant
- US Patent #9119611 | #D783856

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit most InstruSafe Trays

HOW TO ORDER:


EXAMPLE: TRAY IN-2890

Wrapped Foot IN-2890-R	Container Foot IN-2890-TF
---------------------------	------------------------------

FOOT OPTION IS ADDED TO THE END OF THE PART NUMBER

Preconfigured Tray for Arthroscopy

This three-level tray is equipped for both small and large arthroscopic instrument sets. Our tray is sized to fit most three-quarter size containers or can be wrapped for your convenience.



Custom Configured Tray Example for Arthroscopy Example

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES


Cover




Shelf



Basket



Base



ARTHROSCOPY INSTRUMENT TRAY

Arthroscopy Instrument Tray

Holds:
A combination of up to 12 of the following instruments:

- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)	387 L x 244 W x 146 H (mm)
-----------------------------	----------------------------

Wrapped Foot
IN-2846-R

Container Foot
IN-2846-TF

Expand capacity with positionable Hold-Its and Hold-Downs™ silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

CUSTOM CONFIGURED




Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.



CALL US TODAY FOR CUSTOMIZATION!

ARTHROSCOPY TMJ TRAY | 17 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED




AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Tray for Bariatric

The bariatric tray is designed to accommodate even the longest instruments. Our three-level system allows you to get organized and protect all of your bariatric instrumentation in one tray.



Custom Configured Tray Example for Bariatric



BARIATRIC TRAY

Bariatric Instrument Tray

Holds:
A combination of up to 15 (32cm to 45cm+) of the following instruments:


- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars

28 L x 13 W x 8 H (in)	711 L x 330 W x 203 H (mm)
Wrapped Foot IN-8864-R	Container Foot IN-8864-TF

Expand capacity with positionable Hold-Its and Hold-Downs silicone instrument holders.
Multi-level to accommodate a high quantity of instruments.


CALL US TODAY FOR CUSTOMIZATION!

CUSTOM CONFIGURED BARIATRIC TRAY | 31 INSTRUMENTS




BARIATRIC LAPAROSCOPIC SET | 59 INSTRUMENTS

BEFORE




CUSTOM CONFIGURED




AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Cardiovascular















Custom Configured Tray Example for Cardiovascular


MICROVASCULAR SET | 19 INSTRUMENTS

EVANS ROTATABLE TRAY | 32 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED




AFTER

Shelf


Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Dental




Custom Configured Tray Example for Dental


DENTAL SET | 34 INSTRUMENTS

OPERATIVE TRAY SET | 11 INSTRUMENTS

BEFORE




CUSTOM CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

CUSTOM CONFIGURED



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



ENT SYSTEMS



Preconfigured Trays for ENT Modules

Since the mid-1980's, we have offered a complete solution for fine and delicate ear instrumentation through our Middle Ear System. We are proud to claim this as our very first InstruSafe product. Our modules are designed to specifically fit each instrument and can be used individually or in combination, to organize and protect an entire instrument set.

MIX AND MATCH OUR MODULES FROM PAGES 23-24 TO FIT YOUR INSTRUMENT SET, CREATING YOUR UNIQUE ENT SYSTEM!



Tray shown without cover for illustrative purposes.



MIDDLE EAR TRAY

Middle Ear Tray - 2 Level

Holds:

- 2 Alligator Forcep Modules
- 2 Pick + Hook Modules
- 1 Micro-Suction Module
- 1 Speculum Module
- 9½ x 13" Finger Mat

19 5/8 L x 9 5/8 W x 6 H (in)	498 L x 244 W x 152 H (mm)
Wrapped Foot IN-2890-R	Container Foot IN-2890-TF




EAR MODULE TRAY

Alligator Module

Holds:

- 8 Alligator Style Instruments (Cup, Scissors, etc.)

5 7/8 L x 5 1/8 W x 1 H (in)	149 L x 130 W x 25 H (mm)
IN-2001	



EAR MODULE TRAY

Pick + Hook Module, Wide

Holds:

- 10 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 5 1/8 W x 1 H (in)	194 L x 130 W x 25 H (mm)
IN-2002	



EAR MODULE TRAY

Speculum Module

Holds:

- 6-12 Speculas

3 1/2 L x 5 1/8 W x 1 H (in)	89 L x 130 W x 25 H (mm)
IN-2003	

Preconfigured Trays for ENT Modules



Preconfigured Trays for Picks + Hooks



Similar to our original module designs, these standalone trays will protect and organize pick, hook and knife instrumentation for ENT procedures. Instruments can either be positioned at a 45 degree angle for procedural presentation or secured flat for sterilization and storage.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



EAR MODULE TRAY

Micro-Suction Module

Holds:

- 12 Micro-Suction Irrigators with Luer Lock Hubs

5 1/8 L x 3 1/2 W x 1 H (in)	130 L x 89 W x 25 H (mm)
------------------------------	--------------------------

IN-2004



EAR MODULE TRAY

Suction Module

Holds:

- 10 Ferguson Style Suctions

5 7/8 L x 5 1/8 W x 1 H (in)	149 L x 130 W x 25 H (mm)
------------------------------	---------------------------

IN-2005



EAR MODULE TRAY


Alligator Module With Cover

Holds:

- 8 Alligator Style Instruments (Cup, Scissors, etc.)

7 5/8 L x 5 1/8 W x 3 H (in)	194 L x 130 W x 76 H (mm)
------------------------------	---------------------------

IN-2006



EAR MODULE TRAY

Pick + Hook Module

Holds:

- 9 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 4 5/8 W x 1 H (in)	194 L x 117 W x 25 H (mm)
------------------------------	---------------------------

IN-2008

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



PICK + HOOK TRAY


10 Pick + Hook Tray

Holds:

- 10 Pick + Hook Style Instruments

5 1/4 L x 8 3/4 W x 1 1/2 H (in)	133 L x 222 W x 38 H (mm)
----------------------------------	---------------------------

IN-4010



PICK + HOOK TRAY


20 Pick + Hook Tray

Holds:

- 20 Pick + Hook Style Instruments

11 L x 8 3/4 W x 1 1/2 H (in)	279 L x 222 W x 38 H (mm)
-------------------------------	---------------------------

IN-4000



PICK + HOOK TRAY


30 Pick + Hook Tray

Holds:


- 30 Pick + Hook Style Instruments

15 L x 8 3/4 W x 1 1/2 H (in)	381 L x 222 W x 38 H (mm)
-------------------------------	---------------------------

IN-4003



Custom Configured Tray Example for Ear






Custom Configured Tray Example for Ear

MYRINGOTOMY SET | 26 INSTRUMENTS



EAR TRAY | 86 INSTRUMENTS

BEFORE




Shelf

CUSTOM CONFIGURED




AFTER



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





Custom Configured Tray Example for Ear




Preconfigured Tray for Nose/Sinus

MAJOR EAR TRAY | 66 INSTRUMENTS



Shelf



Base




ENDOSCOPY SINUS TRAY

Endoscopy Sinus Instrument Tray

Holds:
A combination of up to 12 of the following instruments:

- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers


15 1/4 L x 9 5/8 W x 5 3/4 H (in)	387 L x 244 W x 146 H (mm)
Wrapped Foot IN-2846-R	Container Foot IN-2846-TF




Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.

CALL US TODAY FOR CUSTOMIZATION!

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Nose/Sinus




Custom Configured Tray Example for Nose/Sinus


SINUS SET | 32 INSTRUMENTS

SINUS TRAY | 34 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED





AFTER



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Shelf

CUSTOM CONFIGURED



Basket




Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Preconfigured Tray for Throat




Custom Configured Tray Example for Throat


Cover




Shelf



Basket



Base



ENDOSCOPY TRAY

Full Size Endoscopy Instrument Tray

Holds:
• 14 Endoscopy Style Instruments

19 7/8 L x 9 5/8 W x 7 H (in)	498 L x 244 W x 178 H (mm)
Wrapped Foot IN-2847-R	Container Foot IN-2847-TF




Small Parts Boxes are perfect for your small delicate instruments or accessories.





Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.

THYROID SET | 25 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED



AFTER

CALL US TODAY FOR CUSTOMIZATION!

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Tray for General

The InstruSafe Battery Tray is designed to accommodate two of the most commonly used sizes of surgical power tool batteries. The tray includes our most trusted medical-grade silicone dividers for the highest level of organization and protection. Our Battery Tray is designed to be wrapped with built-in corner protectors for added security.



BATTERY TRAY

Universal Stryker Battery Tray

Holds:

- 2 Batteries

7 1/2 L x 3 1/2 W x 3 5/8 H (in)	191 L x 89 W x 92 H (mm)
----------------------------------	--------------------------

Wrapped Foot
IN-5015-R

SPECIAL FEATURES


- Silicone divider has inner channels for optimal sterilant penetration
- High durability against automated washing
- Built-In corner protectors for added security
- Lightweight
- FDA 510K Validated




Preconfigured Tray for General

Maximize the protection and organization of your laparoscopic instrument sets with this three-level tray that is excellently constructed to fit both small and large sets of instrumentation.


Cover




Shelf



Basket



Base



GENERAL/GYNECOLOGY TRAY

Extra Long Laparoscopic Instrument Tray

Holds:

A combination of up to 11 of the following instruments:

- | | |
|----------------------------------|------------|
| • Laparoscopic Instruments | • Cannulas |
| • General Dissecting Instruments | • Sheaths |
| • Miscellaneous Parts | • Trocars |
| • Tubing/Cords | |

24 3/4 L x 9 5/8 W x 6 H (in)	629 L x 244 W x 152 H (mm)
-------------------------------	----------------------------

Wrapped Foot
IN-8865-R

Container Foot
IN-8865-TF


Expand capacity with positionable Hold-Its and Hold-Downs silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.




Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.

CALL US TODAY FOR CUSTOMIZATION!



Custom Configured Tray Example for General



Custom Configured Tray Example for General



TENDON STRIPPER | 13 INSTRUMENTS

LAPAROSCOPIC | 47 INSTRUMENTS

BEFORE




CUSTOM CONFIGURED




AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.


Shelf




CUSTOM CONFIGURED




Basket




Base



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for
Hysteroscopy




Custom Configured Tray Example for
Hysteroscopy


HYSTEROSCOPE SET | 10 INSTRUMENTS

HYSTERECTOMY RESECTOSCOPE TRAY | 15 INSTRUMENTS

BEFORE





CUSTOM CONFIGURED




AFTER

CUSTOM CONFIGURED



Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Preconfigured Trays for Kerrisons

Our Neurology and Spine Trays provide maximum accessibility and flexibility when handling and sterilizing kerrisons and forceps.

Preconfigured Trays for Bi-Polar Forceps




Tray shown without cover for illustrative purposes.

KERRISON INSTRUMENT TRAY + RACK

Ten Position Instrument Tray

Holds:
• 10 Kerrisons

15 ¼ L x 9 ⅝ W x 6 H (in)	387 L x 244 W x 152 H (mm)
Wrapped Foot IN-5012-R	Container Foot IN-5012-TF




BAYONET FORCEPS MODULE TRAY

Dual Forceps Module Tray with Cover

Holds:
• 2 Forceps

9 ¾ L x 2 ⅞ W x 1 ⅞ H (in)	248 L x 54 W x 29 H (mm)
IN-5401-02	




KERRISON INSTRUMENT RACK

Ten Position Instrument Rack

Holds:
• 10 Kerrisons

14 ¼ L x 9 ½ W x 6 H (in)	362 L x 241 W x 152 H (mm)
Wrapped Foot N/A	Container Foot IN-5013-TF




BAYONET FORCEPS MODULE TRAY

Dual Extra Long Forceps Module Tray with Cover

Holds:
• 2 Forceps

12 L x 2 ⅞ W x 1 ½ H (in)	305 L x 54 W x 38 H (mm)
IN-5401-03	



KERRISON INSTRUMENT RACK

Five Position Instrument Rack

Holds:
• 5 Kerrisons

14 ¼ L x 5 W x 6 H (in)	362 L x 127 W x 152 H (mm)
Wrapped Foot N/A	Container Foot IN-5014-TF



BAYONET FORCEPS MODULE TRAY

Six Forceps Module Tray with Cover



Holds:
• 6 Forceps

9 ¾ L x 4 ⅝ W x 1 ⅞ H (in)	248 L x 248 W x 29 H (mm)
IN-5401	

Preconfigured Trays for
Bi-Polar Forceps



Custom Configured Tray Example for
Neurology/Spine




BAYONET FORCEPS MODULE TRAY

Eight Forceps Module Tray with Cover

Holds:

- 8 Forceps

12 L x 7 ¼ W x 1 ½ H (in)	305 L x 184 W x 38 H (mm)
IN-5401-08	



BAYONET FORCEPS MODULE TRAY


Twelve Forceps Module Tray with Cover

Holds:


- 12 Forceps

12 L x 9 ⅝ W x 1 ½ H (in)	305 L x 244 W x 38 H (mm)
IN-5401-12	


NEUROLOGY TRAY | 13 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





Custom Configured Tray Example for Neurology/Spine



Custom Configured Tray Example for Neurology/Spine

HIGH SPEED PNEUMATIC DRILL TRAY | 9 INSTRUMENTS

DURAL REPAIR KIT TRAY | 18 INSTRUMENTS




CUSTOM CONFIGURED


CUSTOM CONFIGURED

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Ophthalmology




Custom Configured Tray Example for Ophthalmology


CATARACT-PHACO SET | 28 INSTRUMENTS

EYE TRAY | 57 INSTRUMENTS

BEFORE





CUSTOM CONFIGURED




AFTER

CUSTOM CONFIGURED




Shelf




Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





Custom Configured Tray Example for Ophthalmology



Custom Configured Tray Example for Ophthalmology

CORNEA TRAY | 61 INSTRUMENTS

PHACO HAND PIECE TRAY | 23 INSTRUMENTS




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Preconfigured Trays for Orthopedic

While not usually delicate, orthopedic instrumentation needs to be kept sharp and straight to get the job done. Our Ortho Trays are built to safely accommodate and protect instruments that are commonly found in almost all surgical facilities.



Custom Configured Tray Example for Orthopedic



ORTHOPEDIC TRAY

Lambotte Osteotome Tray

Holds:
 • 12 Instruments—straight or curved, up to 1 1/2" width (38mm)

8 3/8 L x 9 5/8 W x 3 H (in)	213 L x 244 W x 76 H (mm)
Wrapped Foot IN-5010-R	Container Foot IN-5010-TF



CURETTE TRAY


Curette Tray

Holds:
 • 14 Instruments—7-13" in length



15 1/4 L x 9 5/8 W x 2 H (in)	387 L x 244 W x 51 H (mm)
Wrapped Foot IN-5011-R	Container Foot IN-5011-TF

SMALL FRAG SET | 23 INSTRUMENTS

BEFORE




CUSTOM CONFIGURED




AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Plastics



Custom Configured Tray Example for Plastics



MICRO SET | 32 INSTRUMENTS

PLASTIC TRAY | 38 INSTRUMENTS

BEFORE





CUSTOM CONFIGURED



AFTER


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Preconfigured Trays for Bookwalter® Retractors


This tray system is engineered to have a balanced weight distribution that meets the AAMI 25 lb. weight limit standard. We offer a two-level tray design made to accommodate the ring and blade components with an additional tray to store the table post and horizontal bar. Our two-tray system is organized for easy counting and safe for transporting. **Note:** The Mediflex Bookler Retractor also fits in this tray.





Custom Configured Tray Example for Retractors

GREENBERG SET | 17 INSTRUMENTS

BEFORE




CUSTOM CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



BOOKWALTER RETRACTOR TRAY

Ring + Blade Tray

Holds:

- 1 (1") Post Coupling
- 1 Balfour Retractor
- 1 Horizontal Flex Bar
- 1 Harrington Retractor
- 1 Medium Oval Ring
- 1 Segmented Ring
- 1 Gelpi Retractor
- 2 (2x6"/51mm x 152mm) Malleable Retractors
- 1 (3x6"/76mm x 152mm) Malleable Retractor
- 2 (2x3"/51mm x 76mm) Kelly Retractors
- 1 (2x4"/51mm x 110mm) Kelly Retractor
- 1 (2x5"/51mm x 127mm) Kelly Retractor
- 1 (2x6"/51mm x 152mm) Kelly Retractor
- 6 Tilt Ratchet Mechanisms
- 2 Ratchet Mechanisms

22 1/2 L x 13 W x 4 1/2 H (in) 572 L x 330 W x 114 H (mm)

Wrapped Foot
IN-8700-R

Container Foot
IN-8700-TF



BOOKWALTER RETRACTOR TRAY

Post Tray

Holds:


- Table Post and Horizontal Bar

25 L x 5 3/8 W x 3 H (in) 635 L x 137 W x 76 H (mm)


Wrapped Foot
IN-8701-R

Container Foot
IN-8701-TF

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.





Custom Configured Tray Example for
Retractors



Custom Configured Tray Example for
Retractors


BOOKWALTER RETRACTOR TRAY | 38 INSTRUMENTS

CARDIAC MIS RETRACTOR HOLDER TRAY | 2 INSTRUMENTS



Shelf

Base



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Medflex Bookler Retractor is registered trademark of Medflex.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Preconfigured Trays for da Vinci SP®

The storage and reprocessing of the da Vinci SP isn't limited to a single solution. That's why we've designed multiple trays specific to the system's instrumentation to accommodate a wide variety of customer preferences.



Preconfigured Trays for da Vinci SP

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES




INSTRUMENT TRAY

Single Instrument Tray

Holds:

- 1 Instrument

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot IN-8921-R	Container Foot IN-8921-TF




INSTRUMENT TRAY

Dual Instrument Tray

Holds:

- 2 Instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot IN-8922-R	Container Foot IN-8922-TF




INSTRUMENT TRAY

Four Instrument Tray

Holds:

- 4 Instruments

28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8923-R	Container Foot IN-8923-TF




INSTRUMENT TRAY

Six Instrument Tray

Holds:

- 6 Instruments

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8925-R	Container Foot IN-8925-TF




PROCEDURE TRAY

Four Instrument Procedure Tray

Holds:

- 4 Instruments + Associated Accessories

28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8924-R	Container Foot IN-8924-TF



PROCEDURE TRAY

Six Instrument Procedure Tray

Holds:

- 6 Instruments + Associated Accessories

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8926-R	Container Foot IN-8926-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

Preconfigured Trays for da Vinci SP

Preconfigured Trays for da Vinci Xi® Single-Site®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Single-Site instruments—from single instruments to procedural sets and their accommodating accessories—our trays are the recommended solution for instrument reprocessing.



ACCESSORY TRAY

Six Accessory Tray

Holds:

- 2 Cannulas
- 2 Obturators
- 2 Cables

15 ¼ L x 9 ⅝ W x 2 ¾ H (in)	387 L x 244 W x 60 H (mm)
Wrapped Foot IN-8927-R	Container Foot IN-8927-TF



CAMERA TRAY

Single Camera Tray

Holds:

- 1 Camera

28 L x 11 W x 3 H (in)	711 L x 279 W x 76 H (mm)
Wrapped Foot IN-8920-R	Container Foot IN-8920-TF



INSTRUMENT TRAY

Single Instrument Tray

Holds:

- 1 (5mm) Xi Single-Site Instrument

27 ⅜ L x 3 ⅝ W x 3 ½ H (in)	689 L x 86 W x 89 H (mm)
Wrapped Foot IN-8911-R	Container Foot N/A




PROCEDURE TRAY

22 Instrument Procedure Tray

Holds:

- 6 (5mm) Xi Single-Site Instruments
- 3 Rigid Obturators
- 2 Flexible Obturators
- 4 Curved Cannulas
- 3 Rigid Cannulas
- 3 Cables
- 1 Suction Irrigator Tube and Connector

26 ¾ L x 13 W x 3 ⅝ H (in)	679 L x 330 W x 92 H (mm)
Wrapped Foot IN-8912-R	Container Foot N/A




ACCESSORY TRAY

16 Accessory Tray

Holds:

- 3 Rigid Obturators
- 2 Flexible Obturators
- 3 Rigid Cannulas
- 4 Curved Cannulas
- 3 Cables
- 1 Suction Irrigator Tube and Connector

15 ¼ L x 9 ⅝ W x 4 ½ H (in)	387 L x 244 W x 105 H (mm)
Wrapped Foot IN-8913-R	Container Foot IN-8913-TF



All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.




Preconfigured Trays for da Vinci Xi Single-Site



Preconfigured Trays for da Vinci Xi®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Instruments—Single or Dual Staplers and their accommodating accessories—our trays are the recommended solutions for instrument reprocessing.




INSTRUMENT TRAY

Six Instrument Tray

Holds:
• 6 (5mm) Xi Single-Site Instruments

26 3/4 L x 10 W x 3 5/8 H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot IN-8914-R	Container Foot IN-8914-TF

**PAIRS WITH OUR XL CONTAINER
CALL YOUR LOCAL REP TO LEARN MORE!**



INSTRUMENT TRAY

Single Instrument Tray

Holds:
• 1 EndoWrist Instrument

27 1/8 L x 3 1/2 W x 2 1/8 H (in)	689 L x 89 W x 54 H (mm)
Wrapped Foot IN-8931-R	Container Foot N/A




INSTRUMENT TRAY

Two Instrument Tray

Holds:
• 2 EndoWrist instruments

27 1/8 L x 3 3/8 W x 3 1/2 H (in)	689 L x 86 W x 89 H (mm)
Wrapped Foot IN-8947-R	Container Foot N/A



INSTRUMENT TRAY

Three Instrument Tray

Holds:
• 3 EndoWrist instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8948-R	Container Foot IN-8948-TF

All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.


Preconfigured Trays for
da Vinci Xi



Preconfigured Trays for
da Vinci Xi

ARTHRO
BARIATRIC
CARDIO
DENTAL
ENT
GENERAL
GYN
NEURO
EYE
ORTHO
PLASTICS
RETRACTORS
ROBOTICS
SCOPES
TRANSPORT CONTAINERS
ACCESSORIES

ARTHRO
BARIATRIC
CARDIO
DENTAL
ENT
GENERAL
GYN
NEURO
EYE
ORTHO
PLASTICS
RETRACTORS
ROBOTICS
SCOPES
TRANSPORT CONTAINERS
ACCESSORIES



INSTRUMENT TRAY


Four Instrument Tray

Holds:
• 4 EndoWrist Instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in) 679 L x 175 W x 94 H (mm)

Wrapped Foot
IN-8950-R

Container Foot
IN-8950-TF



INSTRUMENT TRAY


Three Instrument + Accessory Tray

Holds:
• 3 EndoWrist Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 3/4 H (in) 679 L x 175 W x 94 H (mm)

Wrapped Foot
IN-8949-R

Container Foot
IN-8949-TF



INSTRUMENT TRAY

Six Instrument Tray


Holds:
• 6 EndoWrist Instruments

26 3/4 L x 10 W x 3 5/8 H (in) 679 L x 254 W x 92 H (mm)

Wrapped Foot
IN-8936-R

Container Foot
IN-8936-TF

**PAIRS WITH OUR XL CONTAINER
CALL YOUR LOCAL REP TO LEARN MORE!**



PROCEDURE TRAY

Four Instrument Procedure Tray

Holds:
• 4 EndoWrist Instruments + Associated Accessories

28 L x 8 W x 5 7/16 H (in) 711 L x 203 W x 138 H (mm)

Wrapped Foot
IN-8915-R

Container Foot
IN-8915-TF

**PAIRS WITH OUR XL CONTAINER
CALL YOUR LOCAL REP TO LEARN MORE!**



INSTRUMENT TRAY


Eight Instrument Tray

Holds:
• 8 EndoWrist Instruments

26 3/4 L x 13 W x 3 5/8 H (in) 679 L x 330 W x 92 H (mm)

Wrapped Foot
IN-8938-R

Container Foot
N/A



PROCEDURE TRAY

Six Instrument Procedure Tray

Holds:
• 6 EndoWrist Instruments + Associated Accessories

28 L x 11 W x 5 3/8 H (in) 711 L x 279 W x 137 H (mm)

Wrapped Foot
IN-8916-R

Tall Foot
IN-8916-TF

**PAIRS WITH OUR XL CONTAINER
CALL YOUR LOCAL REP TO LEARN MORE!**


da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

Preconfigured Trays for
da Vinci Xi



Preconfigured Trays for
da Vinci Xi



PROCEDURE TRAY


Eight Instrument Procedure Tray

Holds:
• 8 EndoWrist Instruments + Associated Accessories

26 3/4 L x 13 W x 3 5/8 H (in) 679 L x 330 W x 92 H (mm)

Wrapped Foot
IN-8935-R

Container Foot
N/A



STAPLER + ACCESSORY TRAY


Single Stapler + Accessory Tray

Holds:
• 1 Stapler Instrument + Associated Accessories

27 1/8 L x 3 3/8 W x 3 1/2 H (in) 689 L x 86 W x 89 H (mm)

Wrapped Foot
IN-8933-R

Container Foot
N/A



STAPLER + ACCESSORY TRAY

Dual Stapler + Accessory Tray

Holds:
• 2 Stapler Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 11/16 H (in) 679 L x 175 W x 94 H (mm)

Wrapped Foot
IN-8932-R

Container Foot
IN-8932-TF

**PAIRS WITH OUR XL CONTAINER
CALL YOUR LOCAL REP TO LEARN MORE!**



ACCESSORY TRAY

Ten Accessory Tray

Holds:
• 6 Cannulas or Obturators
• 1 Release Tool
• 3 Cables

19 5/8 L x 9 5/8 W x 2 H (in) 498 L x 244 W x 51 H (mm)

Wrapped Foot
IN-8939-R

Container Foot
IN-8939-TF



SCOPE TRAY

Single Scope Tray

Holds:
• 1 (8mm) 0° or 30° Scope

24 7/8 L x 13 1/2 W x 4 3/8 H (in) 632 L x 343 W x 111 H (mm)

Wrapped Foot
INS-8930-R

Container Foot
N/A



CAMERA TRAY

Single Camera Tray

Holds:
• 1 Hand Held Camera

15 1/4 L x 9 5/8 W x 3 H (in) 387 L x 244 W x 76 H (mm)

Wrapped Foot
INS-8934-R

Container Foot
INS-8934-TF




**NEED A PLACE TO SOAK?
SEE PAGE 80 FOR DETAILS**

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

Preconfigured Trays for Senhance™ Surgical Robotic Trays


A sophisticated robotic system deserves a solution just as sophisticated to hold its instruments. That's why TransEnterix collaborated with us, the instrument protection tray experts. This line of trays is not only designed and constructed to fit the Senhance Surgical Robotic System instruments and adapters like a glove, but it is also validated specifically for this instrumentation.




Custom Configured Tray Example for Robotic Instrumentation

ROBOTIC LAPAROSCOPIC SET | 9 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED



AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

INSTRUMENT TRAY

Instrument Tray

Holds:

- 8 Instrument Adapters
- 12 Instrument Shafts
- Cables

19 7/8 L x 9 5/8 W x 4 1/4 H (in)	498 L x 244 W x 105 H (mm)
Wrapped Foot IN-8750-R	Container Foot IN-8750-TF

ULTRASONIC INSTRUMENT TRAY

Instrument Tray

Holds:

- 1 Ultrasonic Adapter
- 1 Ultrasonic Transducer

25 L x 3 3/8 W x 3 H (in)	635 L x 137 W x 76 H (mm)
Wrapped Foot IN-8753-R	Container Foot IN-8753-TF

ENDOSCOPE ADAPTER TRAY

Single Instrument Tray

Holds:

- 1 Endoscope Adapter

19 7/8 L x 9 5/8 W x 2 3/8 H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8751-R	Container Foot IN-8751-TF


ENDOSCOPE ADAPTER TRAY

Dual Instrument Tray

Holds:

- 2 Endoscope Adapters

19 7/8 L x 9 5/8 W x 2 3/8 H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8752-R	Container Foot IN-8752-TF



Senhance is a trademark of TransEnterix Surgical, Inc. Please refer to the Senhance™ Surgical Robotic System IFU for instrument reprocessing.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection. da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for Rigid Scopes


Our rigid scope trays are universal and designed to fit scopes of any size for any procedure.



Preconfigured Trays for Rigid Scopes

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



RIGID SCOPE TRAY 4-7"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 4-7" in working length

11 L x 3 ¼ W x 2 ¾ H (in)	279 L x 83 W x 60 H (mm)
Wrapped Foot IN-6011-R	Container Foot N/A




RIGID SCOPE TRAY 9½-21½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 9½-21½" in working length

25 ¼ L x 3 ¼ W x 2 ¾ H (in)	641 L x 83 W x 60 H (mm)
Wrapped Foot IN-6014-R	Container Foot N/A



RIGID SCOPE TRAY 5-13½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 5-13½" in working length

17 ½ L x 3 ¼ W x 2 ¾ H (in)	445 L x 83 W x 60 H (mm)
Wrapped Foot IN-6012-R	Container Foot N/A




RIGID SCOPE TRAY 11-26"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 11-26" in working length

30 ¼ L x 3 ¼ W x 2 ¾ H (in)	765 L x 83 W x 60 H (mm)
Wrapped Foot IN-6015-R	Container Foot N/A



RIGID SCOPE TRAY 7-16½"

Holds:


- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 7-16½" in working length

20 ¼ L x 3 ¼ W x 2 ¾ H (in)	515 L x 83 W x 60 H (mm)
Wrapped Foot IN-6013-R	Container Foot N/A



Preconfigured Trays for Scope + Cable


Our scope + cable trays allow for faster loading before sterilization while also protecting the delicate fiber optics within.



Custom Configured Tray Example for Scope + Cable

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES


- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



SCOPE + CABLE TRAY 5-6"

- Holds:
- Up to 3 (2-10mm) Scopes
 - 1 Cable
 - Accommodates scopes from 5-6" in working length


11 L x 6 3/8 W x 2 3/8 H (in)	279 L x 162 W x 60 H (mm)
Wrapped Foot IN-6031-R	Container Foot N/A



SCOPE + CABLE TRAY 5-13 1/2"

- Holds:
- Up to 3 (2-10mm) Scopes
 - 1 Cable
 - Accommodates scopes from 5-13 1/2" in working length

20 1/2 L x 5 5/8 W x 2 3/8 H (in)	521 L x 143 W x 60 H (mm)
Wrapped Foot IN-6032-R	Container Foot N/A



SCOPE + CABLE TRAY 6 1/2-17"

- Holds:
- Up to 3 (2-10mm) Scopes
 - 1 Cable
 - Accommodates scopes from 6 1/2-17" in working length


25 L x 5 3/8 W x 2 3/8 H (in)	635 L x 137 W x 60 H (mm)
Wrapped Foot IN-6033-R	Container Foot N/A

RESECTOSCOPE | 11 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Preconfigured Trays for Flexible Scopes

Our durable flexible scope trays are sure to protect your delicate and expensive scopes. Our Small Parts Boxes are perfect to keep all your accessories together and protected.




Preconfigured Tray for Flexible Scopes



FLEXIBLE SCOPE TRAY

- Holds:
- 1 Flexible Scope
 - Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6021-R	Container Foot IN-6021-TF



FLEXIBLE SCOPE TRAY

- Holds:
- 1 Flexible Bronchoscope
 - Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6024-R	Container Foot IN-6024-TF



FLEXIBLE SCOPE TRAY

- Holds:
- 1 Flexible Scope with Attached Light Source
 - Accommodates scopes from 12-27" in working length


19 5/8 L x 9 5/8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6022-R	Container Foot IN-6022-TF




FLEXIBLE SCOPE TRAY

- Holds:
- 1 Flexible Video Scope
 - Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6023-R	Container Foot IN-6023-TF



Custom Configured Tray Example for Flexible Scopes



Custom Configured Tray Example for Flexible Scopes

FLEXIBLE URETEROSCOPE | 2 INSTRUMENTS

VIDEO LARYNGOSCOPE | 1 INSTRUMENT



BEFORE





BEFORE



CUSTOM CONFIGURED



CUSTOM CONFIGURED




AFTER


AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





Custom Configured Tray Example for Flexible Scopes



Custom Configured Tray Example for Flexible Scopes

COLONOSCOPY TRAY | 1 INSTRUMENT

FLEXIBLE SCOPE TRAY | 3 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

XL Transport Containers

Transporting da Vinci SP, Xi scopes, and EndoWrist instruments appropriately and effectively is a difficult task. Use the InstruSafe XL Transport Container to get the job done right! This durable container is the perfect size, enabling you to soak and protect your instruments while they travel from the O.R. to sterile processing. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.

XL Transport Container Special Features

- Highly durable material designed to withstand everyday use
- Cart wash safe—maximum temperature of 230°F
- Economical and creates a dedicated space for horizontally soaking robotic instruments
- Made in the USA



IN-PIC-20-D

XL TRANSPORT CONTAINER

XL Transport Container + Optional Drain

Holds:

- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments and accessories

34 L x 17 ¼ W x 6 H (in)
864 L x 438 W x 152 H (mm)

With Drain IN-PIC-20-D	Without Drain IN-PIC-20
---------------------------	----------------------------



IN-PIC-20-D-LT

XL TRANSPORT CONTAINER SYSTEM

XL Transport Container with Lifting Tray + Optional Drain

Holds:

- Lifting Tray
- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments and accessories

Lifting Tray	Transport Container
27 ½ L x 13 ½ W x 3 ⅞ H (in) 699 L x 343 W x 98 H (mm)	34 L x 17 ¼ W x 6 H (in) 864 L x 438 W x 152 H (mm)

With Drain IN-PIC-20-D-LT	Without Drain IN-PIC-20-LT
------------------------------	-------------------------------

LIFTING TRAY

Single Endoscope Configuration

Holds:

- 1 da Vinci Xi or SP Endoscope and accessories

Multiple Instrument Configuration

Holds:

- Up to 8 da Vinci Xi or SP EndoWrist Instruments and accessories

27 ½ L x 13 ½ W x 3 ⅞ H (in)	699 L x 343 W x 98 H (mm)
------------------------------	---------------------------

IN-PIC-34



Single Endoscope Configuration

Multiple Instrument Configuration


Inner gasket and strong silicone latches reduce spills

Dedicated dry erase space on cover



*Top view of lid

Liquid volume feature to accurately measure detergent solutions




Quick drain feature — no need to tip container to drain



LIFTING TRAY

Single Endoscope Configuration

Holds: 1 da Vinci Xi or SP Endoscope



Multiple Instrument Configuration

Holds: up to 8 da Vinci Xi or SP EndoWrist Instruments

Accessory Basket



Large Handles

Loading Identifiers



CONFIGURATION EXAMPLES

Heavy Duty Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. They are similar to our other transport containers but have a more robust construction and feature stainless steel hinges and latches that can be easily undone with just one hand.

Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. Also, you can easily identify content through the transparent amber color. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.



HEAVY DUTY TRANSPORT CONTAINER

<i>Short Container</i>	
20 1/2 L x 12 5/8 W x 3 3/4 H (in)	521 L x 321 W x 95 H (mm)
IN-PIC-15	



HEAVY DUTY TRANSPORT CONTAINER

<i>Tall Container</i>	
20 L x 12 1/4 W x 6 1/8 H (in)	510 L x 310 W x 156 H (mm)
IN-PIC-16	

SPECIAL FEATURES

- Autoclavable
- Made of PEI material (high-temperature plastic)
- High durability against automated washing
- Transparent amber color for easy instrument identification
- Easy locking, removable cover with removable silicone gasket that can be replaced
- Lightweight
- Combine with Instru-Zyme for optimal pre-cleaning



TRANSPORT CONTAINER

<i>Small, Short Container</i>	
13 5/16 L x 15 5/8 W x 4 5/16 H (in)	345 L x 398 W x 110 H (mm)
IN-PIC-14	



TRANSPORT CONTAINER

<i>Large, Short Container</i>	
20 5/8 L x 15 5/8 W x 5 5/16 H (in)	524 L x 397 W x 135 H (mm)
IN-PIC-12	



TRANSPORT CONTAINER

<i>Small, Tall Container</i>	
13 5/16 L x 15 5/8 W x 7 1/2 H (in)	345 L x 398 W x 190 H (mm)
IN-PIC-13	



TRANSPORT CONTAINER

<i>Large, Tall Container</i>	
20 5/8 L x 15 5/8 W x 7 1/2 H (in)	524 L x 397 W x 191 H (mm)
IN-PIC-11	

Silicone Small Parts Boxes

Don't sweat the small stuff. Made of silicone, these boxes are perfect for keeping extra delicate small parts and accessories organized and protected.

Plaque Options

A variety of plaque colors are available to help identify specific sets by procedure or service type. Affix them to our trays and modules using our standard screws. Common identification uses are for bar code placement, tray name labels, or other custom information. Laser engraving of plaques is provided at an additional cost.



SILICONE SMALL PARTS BOX

3x3" Small Parts Box

3 L x 3 W x ¾ H (in)	76 L x 76 W x 19 H (mm)
IN-2652	




SILICONE SMALL PARTS BOX

3x6" Small Parts Box

6 L x 3 W x ¾ H (in)	152 L x 76 W x 19 H (mm)
IN-2653	

- Made entirely of medical grade silicone
- Maximum protection for small parts
- Tethered cover to prevent accidental loss
- Can be attached to trays or placed on posts for easy removal
- Textured internal surface and proper ventilation for efficient sterilant penetration



Small Parts Box is shown within tray. Tray is shown without cover for illustrative purposes only.

Our in-house laser engraving capabilities allow us to personalize trays with custom identification and/or specifications, such as interior instrument descriptions, illustrations, exterior procedure types, specialty, or facility information.


PLAQUES

Plaque Assembly for Trays

1 L x 4 ¾ W (in) 25 L x 111 (mm)	IN-2550-01	SILVER
1 L x 4 ¾ W (in) 25 L x 111 (mm)	IN-2550-05	BLACK
1 L x 4 ¾ W (in) 25 L x 111 (mm)	IN-2550-15	BLACK

Please indicate desired color of plaque using corresponding number below in place of XX using the color options table.

4 ¾ L x 2 W (in) 111 L x 51 W (mm)	IN-2600-XX	●●●
2 L x 2 ½ W (in) 51 L x 56 W (mm)	IN-2300-XX	●●●



COLOR OPTIONS ●●●

Steam Sterilization		Low-Temp Sterilization	
These anodized aluminum plaques for use with steam only.		Plaques have a coating that specially withstands low-temperature sterilization.	
01	SILVER		N/A
02	RED	■	RED 12
03	GREEN	■	GREEN 13
04	BLUE	■	BLUE 14
05	BLACK	■	BLACK 15
07	PURPLE	■	PURPLE 17
08	ORANGE	■	N/A
10	MINT GREEN	■	N/A

Finger Mats


These textured silicone instrument cushion mats are available in a variety of sizes and textures to fit all of your instrument protection needs.

Instru-Zyme®

This powerful enzyme and detergent formulation is designed to assist in the removal of blood, tissue, mucus and other proteinaceous materials. Use this product with our transport containers for a simple solution to water and detergent.


ARTHRO
BARIATRIC
CARDIO
DENTAL
ENT
GENERAL
GYN
NEURO
EYE
ORTHO
PLASTICS
RETRACTORS
ROBOTICS
SCOPES
TRANSPORT CONTAINERS
ACCESSORIES

ARTHRO
BARIATRIC
CARDIO
DENTAL
ENT
GENERAL
GYN
NEURO
EYE
ORTHO
PLASTICS
RETRACTORS
ROBOTICS
SCOPES
TRANSPORT CONTAINERS
ACCESSORIES




FINE FINGER MAT

<i>Full Sheet</i>	
19 L x 13 W (in)	483 L x 330 W (mm)
IN-2616	



FINE FINGER MAT

<i>Half Sheet</i>	
9 ½ L x 13 W (in)	241 L x 330 W (mm)
IN-2617	



FINE FINGER MAT

<i>Quarter Sheet</i>	
4 ¾ L x 13 W (in)	121 L x 330 W (mm)
IN-2618	



INSTRU-ZYME

Properties: • Specific Gravity: 1.200 • pH: 7.8 • Tropical Fragrance		• Color: Light Blue • Form: Liquid Gel	
6 Bottles/Case		IN-1032	

- Will not stain instruments or equipment made of: metal, plastic, glass or rubber
- Non-corrosive
- Biodegradable
- Sustains moisture during transport
- Provides a ready-to-use, neutral pH enzyme solution
- Reduces heavy instrument transport weight as compared to liquid-filled trays
- Minimizes potential exposure due to splashing
- Spray applicator for even coverage of soiled areas
- No rinsing necessary before automated processing

Rounders® Point Guards

Wrap your trays with confidence. Point Guards are engineered to sit securely on almost any tray corner. Prevent holes and tears in your sterilization wrap caused by surgical trays that are heavy, have sharp corners or are improperly handled.



ROUNDERS POINT GUARDS

80/Bag	IN-7330
--------	---------

- Air cushion tray corner protection
- Help avoid compromised sterility in the O.R.
- Can be used in steam, ethylene oxide and low temperature plasma sterilization
- Reusable and economical
- Patent #9119611
- Non-sterile

MADE IN THE USA

Guardian™ Tray

Simple, small instrument sets now have a new home. Our silicone finger mats with sidewalls provide a safe and pliable environment for them to reside in.



GUARDIAN TRAY

<i>Anodized Aluminum Tray with Finger Mat</i>	
10 L x 6 W x 1 ¼ H (in)	254 L x 152 W x 32 H (mm)
IN-1310	



- Accommodate a wide variety of small instrumentation
- Feature built-in and vented silicone finger mats
- Reduce instrument repair and replacement costs
- Allow maximum sterilant penetration
- Prevent moisture accumulation

WISH LIST



WISH LIST





InstruSafe®
Sterile Processing

FOR US + INTERNATIONAL ORDERS:

Summit Medical LLC
An Innovia Medical Company
815 Vikings Parkway
Suite 100
St. Paul, MN 55121
USA

1.888.229.2875
www.innoviamedical.com/instrusafe

FOR UK ORDERS:

DTR Medical Ltd
An Innovia Medical Company
17 Clarion Court, Enterprise Park
Swansea SA6 8RF
UK

+44 (0) 1792 797910
www.dtrmedical.com

InstruSafe Instrument Protection Trays are proudly made in St. Paul, Minnesota, USA.

Summit Medical LLC is an FDA registered facility, ISO 13485:2016 certified and CMDCAS certified.

InstruSafe Trays are 510(k) cleared for multiple sterilization cycles to give our customers more options.

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

da Vinci, SP, S, Xi, Single-Site and EndoWrist are registered trademarks of Intuitive Surgical, Inc. InstruSafe, Instru-Zyme, Hold-Its, Hold-Downs, Rounders Point Guards, and Guardian are trademarks and registered trademarks of Summit Medical.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.

All brands are trademarks or registered trademarks of their respective owners.

©2023 Innovia Medical. All rights reserved. IS-0009 | RevD | JUNE2023