

XL Transport Containers

Transporting da Vinci® SP™, Xi™ scopes, and EndoWrist® instruments appropriately and effectively is a difficult task. The InstruSafe® XL Transport Container is the perfect size to soak and protect your instruments during transport from the O.R. to the SPD.



IN-PIC-20-D



Dry Erase "Soak Start Time" Space

Biohazard Symbol



Water Volume Feature



XL TRANSPORT CONTAINERS

XL Transport Container

Holds:

- da Vinci SP Scopes and EndoWrist Instruments
- da Vinci Xi Scopes and EndoWrist Instruments

34 L x 17 ¼ W x 6 H (in)

864 L x 438 W x 152 H (mm)

IN-PIC-20

XL Transport Container with Drain

Holds:

- da Vinci SP Scopes and EndoWrist Instruments
- da Vinci Xi Scopes and EndoWrist Instruments

34 L x 17 ¼ W x 6 H (in)

864 L x 438 W x 152 H (mm)

IN-PIC-20-D

SPECIAL FEATURES

- Highly durable material designed to withstand everyday use
- Cart wash safe — maximum temperature of 230°F
- Dedicated dry erase "soak start time" space on cover
- Inner gasket and strong silicone latches reduce spills
- Economical and creates a dedicated space for horizontally soaking robotic instruments
- Quick drain feature — no need to tip container to drain
- Built-in biohazard symbol on cover
- Water volume feature to accurately measure detergent solutions
- Made in the USA

*Features on cover emphasized for illustrative purposes only.

da Vinci®, EndoWrist®, Xi™, and SP™ are all trademarks of Intuitive Surgical, Inc.

M-IS-PS-20 | RevB | MARCH2019