



**INTRUSAFE**

An Innovia Medical Brand

**Product Catalog**

*YOUR TRUSTED PARTNER FOR  
SUPERIOR INSTRUMENT PROTECTION*

## About InstruSafe®

InstruSafe® is the leading brand of surgical instrument protection solutions. From the sterile processing department to the operating room and back, we aim to develop and manufacture premium sterilization trays, containers, and reprocessing accessories. Our priority is to improve surgical outcomes, increase long-term savings, and aid patient safety. We create products with superior quality and standards, keeping patients in mind every step of the way.

InstruSafe® is manufactured under Summit Medical LLC, an Innovia Medical® Company. We are your trusted source for specialty surgical and sterile processing essentials that supplies medical devices to over 70 nations worldwide. We work internationally with healthcare organizations and medical device distributors and have established our reputation for using the highest quality of materials and rigorous inspection protocols. Our world-leading brands, comprising over 550 unique medical products, are trusted daily by healthcare professionals globally, elevating patient care delivery and improving clinical outcomes.

## About Innovia Medical®

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

### **ELEVATING CARE. IMPROVING OUTCOMES.**

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

### **INNOVATION LEADERSHIP.**

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

### **UNSURPASSED PRODUCT QUALITY.**

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

### **RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.**

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.



INSTRUSAFE STERILIZATION TRAY ADVANTAGE	
360° OF PROTECTION	4
INDUSTRY DISADVANTAGES	6
FEATURES + BENEFITS	8
ORDER PROCESS: THE INSTRUSAFE WAY	10
STERILIZATION BARRIERS + FOOT OPTIONS	12
TRAYS BY PROCEDURE	
<span style="color: #C00000;">■</span> ARTHROSCOPY	14
<span style="color: #FFA500;">■</span> BARIATRIC	16
<span style="color: #008080;">■</span> CARDIOVASCULAR	18
<span style="color: #000080;">■</span> DENTAL	20
<span style="color: #FF0066;">■</span> ENT	22
<span style="color: #C00000;">■</span> GENERAL	34
<span style="color: #008080;">■</span> GYNECOLOGY	38
<span style="color: #9999CC;">■</span> NEUROLOGY/SPINE	40
<span style="color: #00BFFF;">■</span> OPHTHALMOLOGY	46
<span style="color: #CC0066;">■</span> ORTHOPEDIC	50
<span style="color: #800000;">■</span> PLASTICS	52
<span style="color: #CC99CC;">■</span> RETRACTORS	54
<span style="color: #FFFF00;">■</span> SCOPES	58
ROBOTIC TRAYS	
<span style="color: #008080;">■</span> DAVINCI	68
<span style="color: #000080;">■</span> SENHANCE	78

TRAY ACCESSORIES	
SILICONE SMALL PARTS BOXES	80
PLAQUE OPTIONS	81
FINGER MATS	82
GUARDIAN™ TRAY	82
INSTRU-ZYME®	83
ROUNDERS POINT GUARDS	83
TRANSPORT CONTAINERS	
XL TRANSPORT CONTAINERS	84
HEAVY DUTY TRANSPORT CONTAINERS	86
TRANSPORT CONTAINERS	87
INSTRUSAFE RIGID CONTAINER ADVANTAGE	
FEATURES + BENEFITS	88
RIGID CONTAINERS	
<span style="color: #008080;">■</span> HALF SIZE CONTAINERS	90
<span style="color: #C00000;">■</span> THREE QUARTER SIZE CONTAINERS	92
<span style="color: #FFA500;">■</span> FULL SIZE CONTAINERS	94
<span style="color: #660066;">■</span> XL SIZE CONTAINERS	96
RIGID CONTAINER ACCESSORIES	
LOCKS, FILTERS + DATA CARDS	98
PLAQUE OPTIONS	99
RIGID CONTAINER RESOURCES	
METRIC MEASUREMENT TABLE	100
TRAY CROSS REFERENCE	101
WISH-LIST	102
RULER	106



**SUBMIT PURCHASE ORDERS TO:**

**US + INTERNATIONAL ORDERS**  
**Summit Medical LLC**  
 815 Vikings Parkway, Ste 100  
 St. Paul, MN 55121 USA  
 1-888-229-2875 | +1 651-789-3939  
 www.instrusafe.com

**UK ORDERS**  
**DTR Medical Ltd**  
 17 Clarion Court, Enterprise Park  
 Swansea SA6 8RF UK  
 +44 (0) 1792 797910  
 www.dtrmedical.com

# The InstruSafe® Advantage

## 360° of Protection

The use of complex and delicate instrument sets is a new reality in which sterile processing and the operating room live. Intricacy of surgical instrumentation continues to increase, and so does the cleaning and sterilization process. Now, healthcare facilities face the challenge of reprocessing instruments safely and efficiently while also meeting the demand for fast turnaround. Keeping these sets organized and protected upon arrival in the operating room creates a challenge for all involved. Relieve yourself of these daily obstacles with InstruSafe Trays—an everyday solution for the intricacies of your instrument sets.

We have spent over 30 years understanding the use cycle and improving outcomes. With the use cycle in mind, we have created a product that eliminates common frustrations associated with instrument sets. The products in this catalog are a compilation of our hard work. As you journey through the pages, let us help you solve your everyday instrument problems.

Summit Medical's most recent validations and product advancements are a testament of our dedication to our customers' satisfaction. We wouldn't have continuing success without the healthcare community we serve. Focus, determination and your support will continue our impact on patient care for decades to come.




*Designed to organize and protect from the point of operating room use to decontamination, our solutions benefit everyone who handles our products.*



Designed to ensure optimal sterilization, InstruSafe Trays have undergone rigorous testing before receiving multiple clearances. Cleared for various sterilization cycles, these reliable trays give hospitals and healthcare facilities even more options when it comes to organizing, sterilizing, and protecting their instruments.

For the most up-to-date list of validated cycles, visit [www.instrusafe.com/validations](http://www.instrusafe.com/validations).



POINT CAMERA TO SEE VALIDATIONS



**DECONTAMINATION**

Instruments are individually protected going into the wash cycle and the trays are designed with large perforations, allowing maximum detergent and water flow.

**TRAY PREPARATION**

Each tray is designed specifically for the instruments within it, providing individual locations for each one. This protects them from harming the technician, from damaging themselves and from puncturing or jeopardizing the integrity of the sterile barrier.

**STERILIZATION**

The large amount of perforations in our trays provide superior sterilant penetration. Many of our silicone parts are also designed with minimal contact points and additional ventilation in each slot that match our tray holes to improve sterilant penetration to each individual instrument.

**TRANSPORTATION**

Soft and durable silicone holders in the base and cover secure instruments to prevent shifting and provide 360 degrees of protection. InstruSafe Trays help to avoid unnecessary damage and to present an organized instrument set every time.

**OPERATING ROOM USE**

Each instrument is easily identifiable and placed in a singular slot for surgical presentation and tracking at the end of a procedure.



# Look Familiar?


Unfortunately, this is what most sterile processing departments have come to know. The familiarity and acceptance of these messy, hazardous baskets makes us cringe! Why, you ask? Because the minimal organization and protection of surgical instruments practiced in many facilities presents numerous risks to the individuals that come in contact with them, as well as to the instruments themselves.





# DISADVANTAGES OF MESSY BASKETS

- POTENTIAL TO COMPROMISE STERILE BARRIER BY POKING THROUGH WRAP □
- DANGEROUS TO TRANSPORT FROM THE O.R. TO THE SPD □
- MINIMAL INSTRUMENT PROTECTION AND ORGANIZATION □
- DIFFICULT TO COUNT AND TRACK INSTRUMENTS □
- BASKETS NEED FREQUENT REPLACEMENT □
- ADDS TIME TO SURGICAL SET-UP □
- BROKEN INSTRUMENTS IF DROPPED □
- HIGH RISK FOR FINGER POKES □
- POOR SURGICAL PRESENTATION □
- HARD TO FIND WHAT YOU NEED □
- MISSING/LOST INSTRUMENTS □
- UNNECESSARY REPROCESSING □
- DAMAGED INSTRUMENTS □
- WASTE OF SPACE □




## EXPECT MORE.

For these reasons and many more, InstruSafe Trays offer custom solutions to meet each customer's specific needs. We think it's time for everyone to expect more from their trays.

# Absolute Protection


## Trusted Solutions + Superior Results

Protection starts on the inside. Our comprehensive approach to maximizing instrument longevity is the focal point of our system. We paired our InstruSafe instrument protection trays and our container systems to deliver an unparalleled solution in surgical instrument care.




### QUARTER TURN LATCH

Secures the tray cover to the tray base, providing visual indication of the open/close position. Our off-set latch hole aligns the cover in the correct position ensuring the silicone holders match up with the instruments below.




### RETRACTABLE HANDLES

Our multi-functional space-saving design fits into rigid containers making it effortless to remove the tray.




### TF FOOT

Our Tall Foot fits into rigid containers providing clearance for the filter in the base.



### R FOOT

Rounders Corner Guards prevent tears in sterilization wrap during storage and transportation.










**PERFORATIONS (GRID PATTERN)**

Each tray contains a large quantity of perforations that provide maximum sterilant penetration and ventilation.



**VENTILATION**




Our unique silicone design has channels that align with our perforations allowing detergent and sterilant to flow through the holder and access the entire instrument.



**SILICONE HOLDERS**

InstruSafe silicone holders are soft yet durable, with a history of longevity. Our variety of holders accommodate instruments large or small, allowing each instrument to have its own location within the tray.

Holders in the tray base and cover provide 360 degrees of protection reducing instrument repair and replacement while delivering instruments to the O.R. in their optimal working condition.



**REDUCED CONTACT POINTS**

Each silicone holder is individually designed with ridges and channels built to reduce contact with the instrument while securely holding the instrument in place.

# 5 Simple Steps

If one of our preconfigured tray designs does not fit your needs let our product experts design a custom tray with our effortless process below.




## 1 ASSESS

An InstruSafe Representative conducts a free, on-site assessment of your instrument set.



## 2 GATHER

Your instrument set specifics and inputs are gathered by our team of experts.



## 3 DESIGN

Using your inputs, our product specialists build a tray design with you and your instrument set in mind.







**EVALUATE**

You be the judge! Evaluate the new tray design while it's loaded with your instrument set. Let us know your thoughts and we will make any and all adjustments necessary!



**COMPLETE**

It's all yours! Once complete, each custom tray is assigned its own unique part number, making it simple and easy for you to reorder in future.

**HOW TO USE THE CATALOG**


The colored headers identify which procedure the instrument set you're looking at is used for. The trays on each page are differentiated by preconfigured or custom configured with the symbols shown in the key below. Some trays have several layers including a cover, shelves, baskets, and a base. The custom configured tray examples included in the catalog are a small representation of the trays our local representatives have built for facilities around the world.

**CUSTOM + PRECONFIGURED SYMBOLS**



**THIS SYMBOL IDENTIFIES A PRECONFIGURED TRAY**

Order by part number from your local InstruSafe representative or call us today.



**THIS SYMBOL IDENTIFIES A CUSTOM CONFIGURED TRAY EXAMPLE**

Call us today to schedule your free assessment.



# Sterilization Barriers

InstruSafe Instrument Protection Trays are validated for several sterilization cycles using wrap and/or rigid containers.

For the most up-to-date list of Validated cleared cycles, visit [www.instrusafe.com](http://www.instrusafe.com).

# Foot Options

**R Foot:** Rounders Corner Guards—to be used with wrap

**TF Foot:** Tall Foot fits into rigid containers

STERILE BARRIER	FOOT OPTION
<p><b>WRAP</b></p> 	<p>Rounders Corner Guards <b>-R</b></p> 
<p><b>CONTAINER</b></p> 	<p>Tall Foot <b>-TF</b></p> 

## HOW TO ORDER:

EXAMPLE: TRAY IN-2890	
<p>Wrapped Foot <b>IN-2890-R</b></p>	<p>Container Foot <b>IN-2890-TF</b></p>

FOOT OPTION IS ADDED TO THE END OF THE PART NUMBER

# Rounders<sup>®</sup> Corner Guards

Rounded, medical grade corner guards affix to InstruSafe Tray corners to reduce potential damage to the sterile wrap barrier. Rounders are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to all InstruSafe Trays.



ROUNDERS CORNER GUARDS	
KIT OF FOUR	IN-7334

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit all InstruSafe Trays
- Rounded, medical grade silicone corner guard
- Help avoid compromised sterility in the O.R.
- Environmentally and economically efficient
- Impact resistant
- US Patent #9119611 | #D783856


# Tall Foot

Designed to affix to trays that are intended to be sterilized within a rigid container. Tall Feet are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to most InstruSafe Trays.



TF FOOT KIT	
KIT OF FOUR	IN-7338


- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit most InstruSafe Trays



## Preconfigured Tray for Arthroscopy

This three-level tray is equipped for both small and large arthroscopic instrument sets. Our tray is sized to fit most three-quarter size containers or can be wrapped for your convenience.


Cover




Shelf



Basket



Base



### ARTHROSCOPY TRAY

Holds:  
A combination of up to 12 of the following instruments:

- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)

387 L x 244 W x 146 H (mm)

Wrapped Foot  
IN-2846-R

Container Foot  
IN-2846-TF

Expand capacity with positionable Hold-Its™ and Hold-Downs™ silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

**CUSTOM CONFIGURED**



Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.

**CALL US TODAY FOR CUSTOMIZATION!**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS








## Custom Configured Tray Example for Arthroscopy Example

ARTHROSCOPY TMJ TRAY | 17 INSTRUMENTS

**BEFORE**



**CUSTOM CONFIGURED**



**AFTER**

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



## Preconfigured Tray for Bariatric

The bariatric tray is designed to accommodate even the longest instruments. Our three-level system allows you to get organized and protect all of your bariatric instrumentation in one tray.



### BARIATRIC TRAY

Holds:  
A combination of up to 15 (32cm to 45cm+) of the following instruments:

- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars


28 L x 13 W x 8 H (in)	711 L x 330 W x 203 H (mm)
Wrapped Foot IN-8864-R	Container Foot IN-8864-TF

Expand capacity with positionable Hold-Its™ and Hold-Downs™ silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

## CALL US TODAY FOR CUSTOMIZATION!

### CUSTOM CONFIGURED BARIATRIC TRAY | 31 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS







# Custom Configured Tray Example for Bariatric

BARIATRIC LAPAROSCOPIC SET | 59 INSTRUMENTS

## BEFORE



CUSTOM CONFIGURED



## AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES


ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





## Custom Configured Tray Example for Cardiovascular

ARTHO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES


ROBOTICS

ACCESSORIES



TRANSPORT CONTAINERS

### MICROVASCULAR SET | 19 INSTRUMENTS

# BEFORE




CUSTOM CONFIGURED



# AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.








## Custom Configured Tray Example for Cardiovascular

EVANS ROTATABLE TRAY | 32 INSTRUMENTS

CUSTOM CONFIGURED



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Dental


- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL**
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS

DENTAL SET | 34 INSTRUMENTS

## BEFORE




**CUSTOM CONFIGURED**



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Custom Configured Tray Example for Dental

OPERATIVE TRAY SET | 11 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.


# ENT SYSTEMS

Since the mid-1980's, we have offered a complete solution for fine and delicate ear instrumentation through our Middle Ear System. We are proud to claim this as our very first InstruSafe product. Our modules are designed to specifically fit each instrument and can be used individually or in combination, to organize and protect an entire instrument set.


**MIX AND MATCH OUR MODULES FROM PAGES 23-24 TO FIT YOUR INSTRUMENT SET, CREATING YOUR UNIQUE ENT SYSTEM!**




IN-2004




IN-2001



IN-2003



IN-2008



*Tray shown without cover for illustrative purposes.*

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS





## Preconfigured Trays for ENT Modules



### TWO LEVEL


- Holds:
- 2 Alligator Forcep Modules
  - 2 Pick + Hook Modules
  - 1 Micro-Suction Module
  - 1 Speculum Module
  - 9½ x 13" Finger Mat

19 5/8 L x 9 5/8 W x 6 H (in)

498 L x 244 W x 152 H (mm)

Wrapped Foot  
IN-2890-R

Container Foot  
IN-2890-TF




### ALLIGATOR

- Holds:
- 8 Alligator Style Instruments (Cup, Scissors, etc.)

5 7/8 L x 5 1/8 W x 1 H (in)

149 L x 130 W x 25 H (mm)

IN-2001



### PICK + HOOK

- Holds:
- 10 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 5 1/8 W x 1 H (in)

194 L x 130 W x 25 H (mm)

IN-2002



### SPECULUM

- Holds:
- 6-12 Speculas

3 1/2 L x 5 1/8 W x 1 H (in)

89 L x 130 W x 25 H (mm)

IN-2003

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for ENT Modules

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



## MIRCO SUCTION

Holds:  
• 12 Micro-Suction Irrigators with Luer Lock Hubs

5 1/8 L x 3 1/2 W x 1 H (in)

130 L x 89 W x 25 H (mm)

IN-2004



## SUCTION

Holds:  
• 10 Ferguson Style Suctions

3 1/2 L x 5 1/8 W x 1 H (in)

89 L x 130 W x 25 H (mm)

IN-2005




## ALLIGATOR WITH COVER

Holds:  
• 8 Alligator Style Instruments (Cup, Scissors, etc.)

7 5/8 L x 5 1/8 W x 3 H (in)

194 L x 130 W x 76 H (mm)

IN-2006



## PICK + HOOK

Holds:  
• 9 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 4 5/8 W x 1 H (in)


194 L x 117 W x 25 H (mm)

IN-2008



## Preconfigured Trays for Picks + Hooks

Similar to our original module designs, these standalone trays will protect and organize pick, hook and knife instrumentation for ENT procedures. Instruments can either be positioned at a 45 degree angle for procedural presentation or secured flat for sterilization and storage.




### 10 PICK + HOOK WITH COVER

Holds:  
• 10 Pick + Hook Style Instruments

5 ¼ L x 8 ¾ W x 1 ½ H (in)

133 L x 222 W x 38 H (mm)

IN-4010




### 20 PICK + HOOK WITH COVER

Holds:  
• 20 Pick + Hook Style Instruments

11 L x 8 ¾ W x 1 ½ H (in)

279 L x 222 W x 38 H (mm)

IN-4000



### 30 PICK + HOOK WITH COVER

Holds:  
• 30 Pick + Hook Style Instruments

15 L x 8 ¾ W x 1 ½ H (in)

381 L x 222 W x 38 H (mm)

IN-4003

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES


ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS




Custom Configured Tray Example for  
**Ear**





MYRINGOTOMY SET | 26 INSTRUMENTS

**BEFORE**




**CUSTOM  
CONFIGURED**



**AFTER**

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.







## Custom Configured Tray Example for Ear

EAR TRAY | 86 INSTRUMENTS

CUSTOM CONFIGURED



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS


ACCESSORIES

TRANSPORT CONTAINERS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




# Custom Configured Tray Example for Ear



## MAJOR EAR TRAY | 66 INSTRUMENTS



Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS



## Preconfigured Tray for Nose/Sinus



### ENDOSCOPY SINUS TRAY

Holds:  
A combination of up to 12 of the following instruments:


- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)

387 L x 244 W x 146 H (mm)

Wrapped Foot  
IN-2846-R

Container Foot  
IN-2846-TF



Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.

**CALL US TODAY FOR CUSTOMIZATION!**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS








# Custom Configured Tray Example for Nose/Sinus

SINUS SET | 32 INSTRUMENTS

## BEFORE




CUSTOM CONFIGURED



## AFTER



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT**
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS




## Custom Configured Tray Example for Nose/Sinus


SINUS TRAY | 34 INSTRUMENTS



Basket



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Preconfigured Tray for Throat

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Cover



## FULL SIZE ENDOSCOPY TRAY

Holds:  
• 14 Endoscopy Style Instruments


19 5/8 L x 9 5/8 W x 7 H (in)

498 L x 244 W x 178 H (mm)


Wrapped Foot  
IN-2847-R

Container Foot  
IN-2847-TF


Shelf



Basket



Base




Small Parts Boxes are perfect for your small delicate instruments or accessories.

**CUSTOM CONFIGURED**



Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.


**CALL US TODAY FOR CUSTOMIZATION!**





## Custom Configured Tray Example for Throat

THYROID SET | 25 INSTRUMENTS

**BEFORE**



**CUSTOM CONFIGURED**



**AFTER**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




## Preconfigured Tray for General

Our battery tray options - single-position, three-position, and enhanced two-position - are redesigned to fit various battery sizes seamlessly, ensuring secure organization and protection for surgical power instruments during reprocessing. InstruSafe ensures instrument reliability, enhancing patient safety and procedural outcomes.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL**
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS



UNIVERSAL BATTERY TRAY	
Holds: • 1 Battery	
3 3/4 L x 3 1/2 W x 3 5/8 H (in)	95 L x 89 W x 92 H (mm)
Wrapped Foot IN-5017-R	



UNIVERSAL BATTERY TRAY	
Holds: • 2 Batteries	
7 1/2 L x 3 1/2 W x 3 5/8 H (in)	191 L x 89 W x 92 H (mm)
Wrapped Foot IN-5015-R	



UNIVERSAL BATTERY TRAY	
Holds: • 3 Batteries	
12 L x 3 1/2 W x 3 5/8 H (in)	305 L x 89 W x 92 H (mm)
Wrapped Foot IN-5016-R	







## Preconfigured Tray for General

Maximize the protection and organization of your laparoscopic instrument sets with this three-level tray that is excellently constructed to fit both small and large sets of instrumentation.


Cover




Shelf



Basket



Base



### EXTRA LONG LAPAROSCOPIC TRAY

Holds:  
A combination of up to 11 of the following instruments:

- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars

24 3/4 L x 9 5/8 W x 6 H (in)

629 L x 244 W x 152 H (mm)

Wrapped Foot  
IN-8865-R

Container Foot  
IN-8865-TF

Expand capacity with positionable Hold-Its™ and Hold-Downs™ silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

CUSTOM CONFIGURED



Silicone in basket and base can be customized to accommodate specific instrumentation. Silicone in image is for illustrative purposes only.

**CALL US TODAY FOR CUSTOMIZATION!**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES



TRANSPORT CONTAINERS



# Custom Configured Tray Example for General

TENDON STRIPPER | 13 INSTRUMENTS

## BEFORE




## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL**
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS






# Custom Configured Tray Example for General


LAPAROSCOPIC | 47 INSTRUMENTS



Shelf



Basket



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Custom Configured Tray Example for Hysteroscopy

## HYSTEROSCOPE SET | 10 INSTRUMENTS

# BEFORE



CUSTOM  
CONFIGURED



# AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS







## Custom Configured Tray Example for Hysteroscopy

HYSTERECTOMY RESECTOSCOPE TRAY | 15 INSTRUMENTS

CUSTOM  
CONFIGURED



Shelf

Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



## Preconfigured Trays for Kerrisons

Our Neurology and Spine Trays provide maximum accessibility and flexibility when handling and sterilizing kerrisons and forceps.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



Tray shown without cover for illustrative purposes.

### 10 POSITION INSTRUMENT TRAY + RACK


Holds:  
• 10 Kerrisons

15 ¼ L x 9 ½ W x 6 H (in)

387 L x 244 W x 152 H (mm)

Wrapped Foot  
IN-5012-R

Container Foot  
IN-5012-TF




### 10 POSITION INSTRUMENT RACK

Holds:  
• 10 Kerrisons

14 ¼ L x 9 ½ W x 6 H (in)

362 L x 241 W x 152 H (mm)

Container Foot  
IN-5013-TF



### 5 POSITION INSTRUMENT RACK

Holds:  
• 5 Kerrisons

14 ¼ L x 5 W x 6 H (in)


362 L x 127 W x 152 H (mm)

Container Foot  
IN-5014-TF





## Preconfigured Trays for Bi-Polar Forceps




### BAYONET FORCEPS MODULE TRAY

Holds:  
• 2 Forceps

9 ¾ L x 2 ½ W x 1 ½ H (in)

248 L x 54 W x 29 H (mm)

IN-5401-02




### BAYONET EXTRA LONG FORCEPS MODULE TRAY

Holds:  
• 2 Forceps

12 L x 2 ½ W x 1 ½ H (in)

305 L x 54 W x 38 H (mm)

IN-5401-03



### BAYONET FORCEPS MODULE TRAY

Holds:  
• 6 Forceps

9 ¾ L x 4 ⅝ W x 1 ½ H (in)

248 L x 248 W x 29 H (mm)

IN-5401

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for Bi-Polar Forceps

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS




## BAYONET FORCEPS MODULE TRAY

Holds:  
• 8 Forceps

12 L x 7 ¼ W x 1 ½ H (in)

305 L x 184 W x 38 H (mm)

IN-5401-08



## BAYONET FORCEPS MODULE TRAY


Holds:  
• 12 Forceps

12 L x 9 ⅝ W x 1 ½ H (in)

305 L x 244 W x 38 H (mm)

IN-5401-12







# Custom Configured Tray Example for Neurology/Spine

NEUROLOGY TRAY | 13 INSTRUMENTS

## BEFORE



CUSTOM CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES



TRANSPORT CONTAINERS



## Custom Configured Tray Example for Neurology/Spine


### HIGH SPEED PNEUMATIC DRILL TRAY | 9 INSTRUMENTS

CUSTOM  
CONFIGURED



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.






# Custom Configured Tray Example for Neurology/Spine

## DURAL REPAIR KIT TRAY | 18 INSTRUMENTS

CUSTOM CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




# Custom Configured Tray Example for Ophthalmology

## CATARACT-PHACO SET | 28 INSTRUMENTS

### BEFORE



CUSTOM CONFIGURED




### AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE**
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS







# Custom Configured Tray Example for Ophthalmology

EYE TRAY | 57 INSTRUMENTS

**CUSTOM CONFIGURED**



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES


ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.






# Custom Configured Tray Example for Ophthalmology

## CORNEA TRAY | 61 INSTRUMENTS



Shelf




Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE**
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS





# Custom Configured Tray Example for Ophthalmology

## PHACO HAND PIECE TRAY | 23 INSTRUMENTS



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

YE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



## Preconfigured Trays for Orthopedic

While not usually delicate, orthopedic instrumentation needs to be kept sharp and straight to get the job done. Our Ortho Trays are built to safely accommodate and protect instruments that are commonly found in almost all surgical facilities.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



### LAMBOTTE OSTEOTOME TRAY

Holds:

- 12 Instruments—straight or curved, up to 1 ½" width (38mm)

8 ¾ L x 9 ⅝ W x 3 H (in)	213 L x 244 W x 76 H (mm)
Wrapped Foot IN-5010-R	Container Foot IN-5010-TF




### CURETTE TRAY

Holds:

- 14 Instruments—7-13" in length

15 ¼ L x 9 ⅝ W x 2 H (in)	387 L x 244 W x 51 H (mm)
Wrapped Foot IN-5011-R	Container Foot IN-5011-TF




### KNEE POSITIONER TRAY

Holds:

- Surgical Knee positioning leg holder for knee replacement surgeries.
- 30" Length Base with Carriage, Clamp and Boot

31 ⅞ L x 3 ⅝ W x 6 H (in)	810 L x 270 W x 102 H (mm)
Wrapped Foot IN-5020-R	








## Custom Configured Tray Example for Orthopedic

SMALL FRAG SET | 23 INSTRUMENTS

**BEFORE**



**CUSTOM CONFIGURED**



**AFTER**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS


ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




# Custom Configured Tray Example for Plastics





MICRO SET | 32 INSTRUMENTS

## BEFORE



CUSTOM  
CONFIGURED




## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS**
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS







## Custom Configured Tray Example for Plastics

PLASTIC TRAY | 38 INSTRUMENTS

CUSTOM  
CONFIGURED



ARTHO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Preconfigured Trays for Bookwalter® Retractors

This tray system is engineered to have a balanced weight distribution that meets the AAMI 25 lb. weight limit standard. We offer a two-level tray design made to accommodate the ring and blade components with an additional tray to store the table post and horizontal bar. Our two-tray system is organized for easy counting and safe for transporting. **Note:** The Mediflex Bookler Retractor also fits in this tray.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



## RING + BLADE TRAY

Holds:

- 1 (1") Post Coupling
- 1 Balfour Retractor
- 1 Horizontal Flex Bar
- 1 Harrington Retractor
- 1 Medium Oval Ring
- 1 Segmented Ring
- 1 Gelpi Retractor
- 2 (2x6"/51mm x 152mm) Malleable Retractors
- 1 (3x6"/76mm x 152mm) Malleable Retractor
- 2 (2x3"/51mm x 76mm) Kelly Retractors
- 1 (2x4"/51mm x 110mm) Kelly Retractor
- 1 (2x5"/51mm x 127mm) Kelly Retractor
- 1 (2x6"/51mm x 152mm) Kelly Retractor
- 6 Tilt Ratchet Mechanisms
- 2 Ratchet Mechanisms

22 ½ L x 13 W x 4 ½ H (in)

572 L x 330 W x 114 H (mm)

Wrapped Foot  
IN-8700-R

Container Foot  
IN-8700-TF



## POST TRAY

Holds:

- Table Post and Horizontal Bar

25 L x 5 ⅜ W x 3 H (in)


635 L x 137 W x 76 H (mm)

Wrapped Foot  
IN-8701-R

Container Foot  
IN-8701-TF


The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.






## Custom Configured Tray Example for Retractors

### BOOKWALTER RETRACTOR TRAY | 38 INSTRUMENTS



Shelf



Base


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES

TRANSPORT CONTAINERS





# Custom Configured Tray Example for Retractors

ARTHRO

## GREENBERG SET | 17 INSTRUMENTS

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES



ROBOTICS

ACCESSORIES

# BEFORE




CUSTOM CONFIGURED



# AFTER

TRANSPORT CONTAINERS  
 Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.







## Custom Configured Tray Example for Retractors

### CARDIAC MIS RETRACTOR HOLDER TRAY | 2 INSTRUMENTS

CUSTOM CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



## Preconfigured Trays for Rigid Scopes

Our rigid scope trays are universal and designed to fit scopes of any size for any procedure.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



### RIGID SCOPE TRAY - 4-7"


Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 4-7" in working length

11 L x 3 ¼ W x 2 ⅜ H (in)

279 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6011-R



### RIGID SCOPE TRAY - 5-13½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 5-13½" in working length

17 ½ L x 3 ¼ W x 2 ⅜ H (in)

445 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6012-R



### RIGID SCOPE TRAY - 7-16½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 7-16½" in working length

20 ¼ L x 3 ¼ W x 2 ⅜ H (in)

515 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6013-R

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.





## Preconfigured Trays for Rigid Scopes



### RIGID SCOPE TRAY - 9 1/2-21 1/2"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 9 1/2-21 1/2" in working length

25 1/4 L x 3 1/4 W x 2 3/8 H (in)

641 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6014-R



### RIGID SCOPE TRAY - 11-26"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 11-26" in working length

30 1/8 L x 3 1/4 W x 2 3/8 H (in)

765 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6015-R

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



## Preconfigured Trays for Scope + Cable

Our scope + cable trays allow for faster loading before sterilization while also protecting the delicate fiber optics within.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



### RIGID SCOPE + CABLE TRAY - 5-6"


Holds:

- Up to 3 (2-10mm) Scopes
- 1 Cable
- Accommodates scopes from 5-6" in working length

11 L x 6 3/8 W x 2 3/8 H (in)

279 L x 162 W x 60 H (mm)

Wrapped Foot  
IN-6031-R



### RIGID SCOPE + CABLE TRAY - 5-13 1/2"


Holds:

- Up to 3 (2-10mm) Scopes
- 1 Cable
- Accommodates scopes from 5-13 1/2" in working length

20 1/2 L x 5 5/8 W x 2 3/8 H (in)

521 L x 143 W x 60 H (mm)

Wrapped Foot  
IN-6032-R



### RIGID SCOPE + CABLE TRAY - 6 1/2-17"

Holds:

- Up to 3 (2-10mm) Scopes
- 1 Cable
- Accommodates scopes from 6 1/2-17" in working length


25 L x 5 3/8 W x 2 3/8 H (in)

635 L x 137 W x 60 H (mm)

Wrapped Foot  
IN-6033-R

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.







# Custom Configured Tray Example for Scope + Cable

RESECTOSCOPE | 11 INSTRUMENTS

## BEFORE



CUSTOM CONFIGURED



## AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Resectoscope with three wave resection system and custom instruments for minimally invasive surgery, with a cover to provide 360 degrees of protection.



## Preconfigured Trays for Flexible Scopes

Our durable flexible scope trays are sure to protect your delicate and expensive scopes. Our Small Parts Boxes are perfect to keep all your accessories together and protected.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



### FLEXIBLE SCOPE TRAY

Holds:

- 1 Flexible Scope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot  
IN-6021-R

Container Foot  
IN-6021-TF



### FLEXIBLE SCOPE + LIGHT SOURCE TRAY

Holds:

- 1 Flexible Scope with Attached Light Source
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot  
IN-6022-R

Container Foot  
IN-6022-TF



### FLEXIBLE VIDEO SCOPE TRAY

Holds:

- 1 Flexible Video Scope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot  
IN-6023-R


Container Foot  
IN-6023-TF

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.





## Preconfigured Tray for Flexible Scopes



### FLEXIBLE BRONCHSCOPE TRAY

Holds:

- 1 Flexible Bronchoscope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6024-R	Container Foot IN-6024-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.





## Custom Configured Tray Example for Flexible Scope

FLEXIBLE URETEROSCOPE | 2 INSTRUMENTS

**BEFORE**




**CUSTOM CONFIGURED**



**AFTER**

da Vinci, SP, Xi, EndoWrist and XiGrip are trademarks of Intuitive Surgical, Inc. All other trademarks are the property of their respective owners. Trays come with a cover to provide 360 degrees of protection.





# Custom Configured Tray Example for Flexible Scope

VIDEO LARYNGOSCOPE | 1 INSTRUMENT

## BEFORE



CUSTOM CONFIGURED



## AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

Instrument with EndoWise for single-site laparoscopy is a registered trademark of Safe-Tra Surgical, with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Flexible Scope

## COLONOSCOPY TRAY | 1 INSTRUMENT

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

SCOPES

ROBOTICS


ACCESSORIES

TRANSPORT CONTAINERS



Tray shown without covers for illustrative purposes. All custom configured trays are marked with a red dot, provide 360 degrees of protection.






## Custom Configured Tray Example for Flexible Scope

### FLEXIBLE SCOPE TRAY | 3 INSTRUMENTS

**CUSTOM CONFIGURED**



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

*Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection. da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.*



## Preconfigured Trays for da Vinci SP®

The storage and reprocessing of the da Vinci SP isn't limited to a single solution. That's why we've designed multiple trays specific to the system's instrumentation to accommodate a wide variety of customer preferences.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES


TRANSPORT CONTAINERS



### SINGLE INSTRUMENT TRAY

Holds:  
• 1 Instrument


26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot IN-8921-R	Container Foot IN-8921-TF



### DUAL INSTRUMENT TRAY

Holds:  
• 2 Instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot IN-8922-R	Container Foot IN-8922-TF



### FOUR INSTRUMENT TRAY


Holds:  
• 4 Instruments

28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8923-R	Container Foot IN-8923-TF






## Preconfigured Trays for da Vinci SP



### SIX INSTRUMENT TRAY

Holds:  
• 6 Instruments


28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8925-R	Container Foot IN-8925-TF



### FOUR INSTRUMENT PROCEDURE TRAY

Holds:  
• 4 Instruments + Associated Accessories

28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8924-R	Container Foot IN-8924-TF



### SIX INSTRUMENT PROCEDURE TRAY

Holds:  
• 6 Instruments + Associated Accessories

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8926-R	Container Foot IN-8926-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for da Vinci SP



## SIX ACCESSORY TRAY

Holds:

- 2 Cannulas
- 2 Obturators
- 2 Cables

15 ¼ L x 9 ⅝ W x 2 ⅜ H (in)	387 L x 244 W x 60 H (mm)
Wrapped Foot IN-8927-R	Container Foot IN-8927-TF



## SINGLE CAMERA TRAY

Holds:

- 1 Camera

28 L x 11 W x 3 H (in)	711 L x 279 W x 76 H (mm)
Wrapped Foot IN-8920-R	Container Foot IN-8920-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS





## Preconfigured Trays for da Vinci Xi® Single-Site®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Single-Site instruments—from single instruments to procedural sets and their accommodating accessories—our trays are the recommended solution for instrument reprocessing.




### Xi SINGLE INSTRUMENT TRAY

- Holds:
- 1 (5mm) Xi Single-Site Instrument

27 1/8 L x 3 3/8 W x 3 1/2 H (in)	689 L x 86 W x 89 H (mm)
-----------------------------------	--------------------------

Wrapped Foot  
IN-8911-R




### Xi 22 INSTRUMENT PROCEDURE TRAY

- Holds:
- 6 (5mm) Xi Single-Site Instruments
  - 3 Rigid Obturators
  - 2 Flexible Obturators
  - 4 Curved Cannulas
  - 3 Rigid Cannulas
  - 3 Cables
  - 1 Suction Irrigator Tube and Connector

26 3/4 L x 13 W x 3 5/8 H (in)	679 L x 330 W x 92 H (mm)
--------------------------------	---------------------------

Wrapped Foot  
IN-8912-R



### Xi 16 ACCESSORY TRAY

- Holds:
- 3 Rigid Obturators
  - 2 Flexible Obturators
  - 3 Rigid Cannulas
  - 4 Curved Cannulas
  - 3 Cables
  - 1 Suction Irrigator Tube and Connector

15 1/4 L x 9 5/8 W x 4 1/8 H (in)	387 L x 244 W x 105 H (mm)
-----------------------------------	----------------------------

Wrapped Foot  
IN-8913-R

Container Foot  
IN-8913-TF

All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES


ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for da Vinci Xi Single-Site



## Xi SIX INSTRUMENT TRAY

Holds:

- 6 (5mm) Xi Single-Site Instruments

26 ¾ L x 10 W x 3 ⅝ H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot IN-8914-R	Container Foot IN-8914-TF

**Pairs with our XL Container**  
**Call your local Rep to learn more!**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

 All da Vinci Xi Single-Site Trays have green covers.

*Excludes IN-8911*





## Preconfigured Trays for da Vinci Xi®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Instruments—Single or Dual Staplers and their accommodating accessories—our trays are the recommended solutions for instrument reprocessing.




### Xi SINGLE INSTRUMENT TRAY

Holds:  
• 1 EndoWrist Instrument

27 1/8 L x 3 1/2 W x 2 1/8 H (in)	689 L x 89 W x 54 H (mm)
-----------------------------------	--------------------------

Wrapped Foot  
IN-8931-R




### Xi DUAL INSTRUMENT TRAY

Holds:  
• 2 EndoWrist instruments

27 1/8 L x 3 3/8 W x 3 1/2 H (in)	689 L x 86 W x 89 H (mm)
-----------------------------------	--------------------------

Wrapped Foot  
IN-8947-R



### Xi THREE INSTRUMENT TRAY

Holds:  
• 3 EndoWrist instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
-----------------------------------	---------------------------

Wrapped Foot  
IN-8948-R

Container Foot  
IN-8948-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for da Vinci Xi

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES


TRANSPORT CONTAINERS



## Xi FOUR INSTRUMENT TRAY

Holds:  
• 4 EndoWrist Instruments

26 ¾ L x 6 ⅞ W x 3 ¾ H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8950-R	Container Foot IN-8950-TF



## Xi SIX INSTRUMENT TRAY

Holds:  
• 6 EndoWrist Instruments

26 ¾ L x 10 W x 3 ⅞ H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot IN-8936-R	Container Foot IN-8936-TF

Pairs with our XL Container  
Call your local Rep to learn more!



## Xi EIGHT INSTRUMENT TRAY


Holds:  
• 8 EndoWrist Instruments

26 ¾ L x 13 W x 3 ⅞ H (in)	679 L x 330 W x 92 H (mm)
Wrapped Foot IN-8938-R	





## Preconfigured Trays for da Vinci Xi




### Xi THREE INSTRUMENT + ACCESSORY TRAY

Holds:  
 • 3 EndoWrist Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
-----------------------------------	---------------------------

Wrapped Foot  
 IN-8949-R

Container Foot  
 IN-8949-TF



### Xi FOUR INSTRUMENT PROCEDURE TRAY


Holds:  
 • 4 EndoWrist Instruments + Associated Accessories

28 L x 8 W x 5 7/16 H (in)	711 L x 203 W x 138 H (mm)
----------------------------	----------------------------

Wrapped Foot  
 IN-8915-R

Container Foot  
 IN-8915-TF

**Pairs with our XL Container**  
**Call your local Rep to learn more!**



### Xi SIX INSTRUMENT PROCEDURE TRAY

Holds:  
 • 6 EndoWrist Instruments + Associated Accessories

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
----------------------------	----------------------------

Wrapped Foot  
 IN-8916-R

Container Foot  
 IN-8916-TF

**Pairs with our XL Container**  
**Call your local Rep to learn more!**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for da Vinci Xi

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



## Xi EIGHT INSTRUMENT PROCEDURE TRAY


Holds:

- 8 EndoWrist Instruments + Associated Accessories

26 3/4 L x 13 W x 3 5/8 H (in)

679 L x 330 W x 92 H (mm)

Wrapped Foot  
IN-8935-R



## SINGLE STAPLER + ACCESSORY TRAY


Holds:

- 1 Stapler Instrument + Associated Accessories

27 1/8 L x 3 3/8 W x 3 1/2 H (in)

689 L x 86 W x 89 H (mm)

Wrapped Foot  
IN-8933-R



## DUAL STAPLER + ACCESSORY TRAY

Holds:

- 2 Stapler Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 1/16 H (in)

679 L x 175 W x 94 H (mm)

Wrapped Foot  
IN-8932-R

Container Foot  
IN-8932-TF

**Pairs with our XL Container  
Call your local Rep to learn more!**





## Preconfigured Trays for da Vinci Xi



### TEN ACCESSORY TRAY


- Holds:
- 6 Cannulas or Obturators
  - 1 Release Tool
  - 3 Cables

19 5/8 L x 9 5/8 W x 2 H (in)

498 L x 244 W x 51 H (mm)

Wrapped Foot  
IN-8939-R

Container Foot  
IN-8939-TF



### SINGLE SCOPE TRAY

- Holds:
- 1 (8mm) 0° or 30° Scope

24 7/8 L x 13 1/2 W x 4 3/8 H (in)

632 L x 343 W x 111 H (mm)

Wrapped Foot  
INS-8930-R



### SINGLE CAMERA TRAY

- Holds:
- 1 Hand Held Camera

15 1/4 L x 9 5/8 W x 3 H (in)

387 L x 244 W x 76 H (mm)

Wrapped Foot  
INS-8934-R

Container Foot  
INS-8934-TF



NEED A PLACE TO SOAK?  
SEE PAGE 84 FOR DETAILS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Preconfigured Trays for *Senhance*<sup>™</sup> Surgical Robotic Trays

A sophisticated robotic system deserves a solution just as sophisticated to hold its instruments. That's why TransEnterix collaborated with us, the instrument protection tray experts. This line of trays is not only designed and constructed to fit the Senhance Surgical Robotic System instruments and adapters like a glove, but it is also validated specifically for this instrumentation.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS





RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



## INSTRUMENT TRAY

Holds:

- 8 Instrument Adapters
- 12 Instrument Shafts
- Cables

19 5/8 L x 9 5/8 W x 4 1/8 H (in)

498 L x 244 W x 105 H (mm)

Wrapped Foot  
IN-8750-R

Container Foot  
IN-8750-TF

## ULTRASONIC INSTRUMENT TRAY

Holds:

- 1 Ultrasonic Adapter
- 1 Ultrasonic Transducer

25 L x 3 3/8 W x 3 H (in)

635 L x 137 W x 76 H (mm)

Wrapped Foot  
IN-8753-R

Container Foot  
IN-8753-TF

## ENDOSCOPE ADAPTER TRAY

Holds:

- 1 Endoscope Adapter

19 5/8 L x 9 5/8 W x 2 3/8 H (in)

498 L x 244 W x 60 H (mm)

Wrapped Foot  
IN-8751-R

Container Foot  
IN-8751-TF

## DUAL ENDOSCOPE ADAPTER TRAY

Holds:

- 2 Endoscope Adapters

19 5/8 L x 9 5/8 W x 2 3/8 H (in)


498 L x 244 W x 60 H (mm)

Wrapped Foot  
IN-8752-R

Container Foot  
IN-8752-TF

*Senhance is a trademark of TransEnterix Surgical, Inc. Please refer to the Senhance<sup>™</sup> Surgical Robotic System IFU for instrument reprocessing.*







# Custom Configured Tray Example for Robotic Instrumentation

ROBOTIC LAPAROSCOPIC SET | 9 INSTRUMENTS

**BEFORE**



**CUSTOM CONFIGURED**



**AFTER**

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

# Silicone Small Parts Boxes

Don't sweat the small stuff. Made of silicone, these boxes are perfect for keeping extra delicate small parts and accessories organized and protected.



## SILICONE SMALL PARTS BOX


3"x3" Small Parts Box	
3 L x 3 W x ¾ H (in)	76 L x 76 W x 19 H (mm)
IN-2652	



## SILICONE SMALL PARTS BOX

3"x6" Small Parts Box	
6 L x 3 W x ¾ H (in)	152 L x 76 W x 19 H (mm)
IN-2653	

- Made entirely of medical grade silicone
- Maximum protection for small parts
- Tethered cover to prevent accidental loss
- Can be attached to trays or placed on posts for easy removal
- Textured internal surface and proper ventilation for efficient sterilant penetration



*Small Parts Box is shown within tray. Tray is shown without cover for illustrative purposes only.*




# Plaque Options





A variety of plaque colors are available to help identify specific sets by procedure or service type. Affix them to our trays and modules using our standard screws. Common identification uses are for bar code placement, tray name labels, or other custom information. Laser engraving of plaques is provided at an additional cost.

Our in-house laser engraving capabilities allow us to personalize trays with custom identification and/or specifications, such as interior instrument descriptions, illustrations, exterior procedure types, specialty, or facility information.

## PLAQUE ASSEMBLY FOR TRAYS

<b>SILVER</b>	IN-2550-01		
	1 L x 4 3/8 W (in)		25 L x 111 (mm)
<b>BLACK</b>	IN-2550-05		
	1 L x 4 3/8 W (in)		25 L x 111 (mm)
<b>BLACK</b>	IN-2550-15		
	1 L x 4 3/8 W (in)		25 L x 111 (mm)
	IN-2600-XX	Please indicate desired color of plaque using corresponding number below in place of XX using the color options table.	
	4 3/8 L x 2 W (in)		111 L x 51 W (mm)
	IN-2300-XX	Please indicate desired color of plaque using corresponding number below in place of XX using the color options table.	
	2 L x 2 1/8 W (in)		51 L x 56 W (mm)



		COLOR OPTIONS		
<b>Steam Sterilization</b>		<b>Low-Temp Sterilization</b>		
These anodized aluminum plaques for use with steam only.		Plaques have a coating that specially withstands low-temperature sterilization.		
01	SILVER		N/A	
02	RED		RED	12
03	GREEN		GREEN	13
04	BLUE		BLUE	14
05	BLACK		BLACK	15
07	PURPLE		PURPLE	17
08	ORANGE		N/A	
10	MINT GREEN		N/A	


ARTHRO  
BARIATRIC  
CARDIO  
DENTAL  
ENT  
GENERAL  
GYN  
NEURO  
EYE  
ORTHO  
PLASTICS  
RETRACTORS  
SCOPES  
ROBOTICS  
ACCESSORIES

TRANSPORT CONTAINERS


# Finger Mats

These textured silicone instrument cushion mats are available in a variety of sizes and textures to fit all of your instrument protection needs.


- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS



FINE FINGER MAT	
Full Sheet	
19 L x 13 W (in)	483 L x 330 W (mm)
IN-2616	



FINE FINGER MAT	
Half Sheet	
9 ½ L x 13 W (in)	241 L x 330 W (mm)
IN-2617	



FINE FINGER MAT	
Quarter Sheet	
4 ¾ L x 13 W (in)	121 L x 330 W (mm)
IN-2618	

## Guardian™ Tray

Simple, small instrument sets now have a new home. Our silicone finger mats with sidewalls provide a safe and pliable environment for them to reside in.



GUARDIAN TRAY	
Anodized Aluminum Tray with Finger Mat	
10 L x 6 W x 1 ¼ H (in)	254 L x 152 W x 32 H (mm)
IN-1310	

- Accommodate a wide variety of small instrumentation
- Feature built-in and vented silicone finger mats
- Reduce instrument repair and replacement costs
- Allow maximum sterilant penetration
- Prevent moisture accumulation



# Instru-Zyme®

This powerful enzyme and detergent formulation is designed to assist in the removal of blood, tissue, mucus and other proteinaceous materials. Use this product with our transport containers for a simple solution to water and detergent.



## FINE FINGER MAT

Properties:

- Specific Gravity: 1.200
- pH: 7.8
- Tropical Fragrance
- Color: Light Blue
- Form: Liquid Gel

Currently for U.S Sales Only


Quantity  
6 Bottles/Case

IN-1032

- Will not stain instruments or equipment made of: metal, plastic, glass or rubber
- Non-corrosive
- Biodegradable
- Sustains moisture during transport
- Provides a ready-to-use, neutral pH enzyme solution
- Reduces heavy instrument transport weight as compared to liquid-filled trays
- Minimizes potential exposure due to splashing
- Spray applicator for even coverage of soiled areas
- No rinsing necessary before automated processing

# Rounders® Point Guards

Wrap your trays with confidence. Point Guards are engineered to sit securely on almost any tray corner. Prevent holes and tears in your sterilization wrap caused by surgical trays that are heavy, have sharp corners or are improperly handled.



## ROUNDERS POINT GUARDS

Reusable and economical

Quantity  
80/Bag

IN-7330

- Air cushion tray corner protection
- Help avoid compromised sterility in the O.R.
- Can be used in steam, ethylene oxide and low temperature plasma sterilization
- Patent #9119611
- Non-sterile

**MADE IN THE USA**

# XL Transport Containers

Transporting da Vinci SP, Xi scopes, and EndoWrist instruments appropriately and effectively is a difficult task. Use the InstruSafe XL Transport Container to get the job done right! This durable container is the perfect size, enabling you to soak and protect your instruments while they travel from the O.R. to sterile processing. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.



- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- SCOPES
- ROBOTICS
- ACCESSORIES
- TRANSPORT CONTAINERS

## XL TRANSPORT CONTAINER

Holds:

- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments and accessories

34 L x 17 ¼ W x 6 H (in)	864 L x 438 W x 152 H (mm)
With Drain IN-PIC-20-D	Without Drain IN-PIC-20

## XL TRANSPORT CONTAINER WITH LIFTING TRAY

Holds:

- Lifting Tray
- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments and accessories

34 L x 17 ¼ W x 6 H (in)	864 L x 438 W x 152 H (mm)
*LIFTING TRAY DIMENSIONS BELOW	
With Drain IN-PIC-20-D-LT	Without Drain IN-PIC-20-LT

## SINGLE ENDOSCOPE LIFTING TRAY

Holds:

- 1 da Vinci Xi or SP Endoscope and accessories

27 ½ L x 13 ½ W x 3 ⅞ H (in)	699 L x 343 W x 98 H (mm)
IN-PIC-34	

## MULTIPLE INSTRUMENT LIFTING TRAY

Holds:

- Up to 8 da Vinci Xi or SP EndoWrist Instruments and accessories

27 ½ L x 13 ½ W x 3 ⅞ H (in)	699 L x 343 W x 98 H (mm)
IN-PIC-34	



# XL Transport Container Special Features

- Highly durable material designed to withstand everyday use
- Cart wash safe—maximum temperature of 230°F
- Economical and creates a dedicated space for horizontally soaking robotic instruments
- Made in the USA

Inner gasket and strong silicone latches reduce spills

Dedicated dry erase space on cover



\*Top view of lid



Liquid volume feature to accurately measure detergent solutions




Quick drain feature — no need to tip container to drain



## LIFTING TRAY

### Single Endoscope Configuration

Holds: 1 da Vinci Xi or SP Endoscope



### Multiple Instrument Configuration

Holds: up to 8 da Vinci Xi or SP EndoWrist Instruments

### Accessory Basket



Large Handles

Loading Identifiers



CONFIGURATION EXAMPLES

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

TRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

# Heavy Duty Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. They are similar to our other transport containers but have a more robust construction and feature stainless steel hinges and latches that can be easily undone with just one hand.



HEAVY DUTY TRANSPORT CONTAINER	
Short Container	
20 1/2 L x 12 5/8 W x 3 3/4 H (in)	521 L x 321 W x 95 H (mm)
IN-PIC-15	



HEAVY DUTY TRANSPORT CONTAINER	
Tall Container	
20 L x 12 1/4 W x 6 1/8 H (in)	510 L x 310 W x 156 H (mm)
IN-PIC-16	

## SPECIAL FEATURES

- Autoclavable
- Made of PEI material (high-temperature plastic)
- High durability against automated washing
- Transparent amber color for easy instrument identification
- Easy locking, removable cover with removable silicone gasket that can be replaced
- Lightweight
- Combine with Instru-Zyme for optimal pre-cleaning: see page 83

ARTHO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES

ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS



# Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. Also, you can easily identify content through the transparent amber color. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.



TRANSPORT CONTAINER	
Small, Short Container	
13 5/16 L x 15 5/8 W x 4 5/16 H (in)	345 L x 398 W x 110 H (mm)
IN-PIC-14	



TRANSPORT CONTAINER	
Large, Short Container	
20 5/8 L x 15 5/8 W x 5 5/16 H (in)	524 L x 397 W x 135 H (mm)
IN-PIC-12	



TRANSPORT CONTAINER	
Small, Tall Container	
13 5/16 L x 15 5/8 W x 7 1/2 H (in)	345 L x 398 W x 190 H (mm)
IN-PIC-13	



TRANSPORT CONTAINER	
Large, Tall Container	
20 5/8 L x 15 5/8 W x 7 1/2 H (in)	524 L x 397 W x 191 H (mm)
IN-PIC-11	

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

SCOPES


ROBOTICS

ACCESSORIES

TRANSPORT CONTAINERS

# Rigid Container Feature + Benefits

We provide sterile processing and the operating room a system that protects investments inside and out. Our unique instrument protection trays and durable external containers are designed to dramatically reduce instrument repairs and replacement costs. This combination delivers uncompromised surgical instrument sets to aid in a higher quality of clinical care and helps maximize economic value through improved instrument set longevity.







### ANODIZED ALUMINUM

Durable, lightweight material, substantially lighter than competitive containers



### FILTER RETENTION PLATES

Easy to install with audible snap when locked in position



### STAINLESS STEEL LATCHES

One step stainless steel locking mechanism eliminates malfunctions



### ERGONOMIC HANDLES

90° stop handles feature ergonomic silicone and flip to read plaques and data cards



### STANDARD + UNIQUE SIZES

Half, three quarters, full and unique XL sizes available



# Lid + Solid Bottoms

HALF

THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES



HALF SIZE LID	
Outside 11 3/4 L x 11 W x 1 1/4 H (in)	
CS-1001	



HALF SIZE 4" SOLID BOTTOM	
Outside with Lid 11 3/4 L x 11 W x 4 1/8 H (in)	Usable Inside 10 1/4 L x 10 1/4 W x 3 1/4 H (in)
CS-1002	



HALF SIZE 6" SOLID BOTTOM	
Outside with Lid 11 3/4 L x 11 W x 5 5/8 H (in)	Usable Inside 10 1/4 L x 10 1/4 W x 4 1/2 H (in)
CS-1003	



# Perforated Bottoms + Basket




## HALF SIZE 4" PERFORATED BOTTOM

Outside with Lid	Usable Inside
11 3/4 L x 11 W x 4 1/8 H (in)	10 1/4 L x 10 1/4 W x 2 3/4 H (in)
CS-1012	



## HALF SIZE 6" PERFORATED BOTTOM

Outside with Lid	Usable Inside
11 3/4 L x 11 W x 5 5/8 H (in)	10 1/4 L x 10 1/4 W x 4 H (in)
CS-1013	



## HALF SIZE CONTAINER BASKET

Outside with Lid	Usable Inside
8 1/2 L x 9 3/4 W x 2 3/4 H (in)	8 3/8 L x 9 5/8 W x 2 H (in)
CS-1120	

HALF SIZE

THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES

# Lid + Solid Bottoms

HALF


THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES



## THREE QUARTER SIZE LID

Outside  
18 ¼ L x 11 W x 1 ¼ H (in)

CS-2001



## THREE QUARTER SIZE 4" SOLID BOTTOM

Outside with Lid  
18 ¼ L x 11 W x 4 ½ H (in)

Usable Inside  
16 ¾ L x 10 ¼ W x 3 ¼ H (in)

CS-2002



## THREE QUARTER SIZE 6" SOLID BOTTOM

Outside with Lid  
18 ¼ L x 11 W x 5 ½ H (in)

Usable Inside  
16 ¾ L x 10 ¼ W x 4 ½ H (in)

CS-2003



# Perforated Bottoms + Basket



## THREE QUARTER SIZE 4" PERFORATED BOTTOM

Outside with Lid  
18 ¼ L x 11 W x 4 ½ H (in)

Usable Inside  
16 ¾ L x 10 ¼ W x 2 ¾ H (in)

CS-2012




## THREE QUARTER SIZE 6" PERFORATED BOTTOM

Outside with Lid  
18 ¼ L x 11 W x 5 ½ H (in)

Usable Inside  
16 ¾ L x 10 ¼ W x 4 H (in)

CS-2013



## THREE QUARTER SIZE CONTAINER BOTTOM

Outside with Lid  
15 ¾ L x 9 ¾ W x 2 ¾ H (in)

Usable Inside  
15 ¼ L x 9 ¾ W x 2 H (in)

CS-2120

HALF SIZE

THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES

# Lid + Solid Bottoms

HALF


THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES



## FULL SIZE LID

Outside 23 1/8 L x 11 W x 1 1/4 H (in)
CS-3001



## FULL SIZE 4" SOLID BOTTOM

Outside with Lid 23 1/8 L x 11 W x 4 1/8 H (in)	Usable Inside 21 5/8 L x 10 1/4 W x 3 1/4 H (in)
CS-3002	



## FULL SIZE 6" SOLID BOTTOM

Outside with Lid 23 1/8 L x 11 W x 5 5/8 H (in)	Usable Inside 21 5/8 L x 10 1/4 W x 4 1/2 H (in)
CS-3003	



## FULL SIZE 8" SOLID BOTTOM

Outside with Lid 23 1/8 L x 11 W x 8 1/4 H (in)	Usable Inside 21 5/8 L x 10 1/4 W x 7 1/8 H (in)
CS-3004	



## FULL SIZE 10" SOLID BOTTOM

Outside with Lid 23 1/8 L x 11 W x 10 3/8 H (in)	Usable Inside 21 5/8 L x 10 1/4 W x 9 1/2 H (in)
CS-3005	

# Perforated Bottoms + Baskets



## FULL SIZE 4" PERFORATED BOTTOM

Outside with Lid  
23 1/8 L x 11 W x 4 1/8 H (in)

Usable Inside  
21 5/8 L x 10 1/4 W x 2 3/4 H (in)

CS-3012



## FULL SIZE 6" PERFORATED BOTTOM

Outside with Lid  
23 1/8 L x 11 W x 5 5/8 H (in)

Usable Inside  
21 5/8 L x 10 1/4 W x 4 H (in)

CS-3013



## FULL SIZE 8" PERFORATED BOTTOM

Outside with Lid  
23 1/8 L x 11 W x 8 1/4 H (in)

Usable Inside  
21 5/8 L x 10 1/4 W x 6 5/8 H (in)

CS-3014



## FULL SIZE 10" PERFORATED BOTTOM

Outside with Lid  
23 1/8 L x 11 W x 10 3/8 H (in)

Usable Inside  
21 5/8 L x 10 1/4 W x 9 H (in)

CS-3015



## FULL SIZE CONTAINER BASKETS

Outside with Lid  
19 3/4 L x 9 3/4 W x 2 3/4 H (in)

Usable Inside  
19 5/8 L x 9 5/8 W x 2 H (in)

CS-3120

Outside with Lid  
19 3/4 L x 9 3/4 W x 6 3/4 H (in)

Usable Inside  
19 5/8 L x 9 5/8 W x 6 H (in)

CS-3160

HALF SIZE

THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES



# Lid, Bottom + Basket

HALF

THREE QUARTER

FULL

XL

ACCESSORIES


RESOURCES



XL SIZE LID	
Outside 30 3/8 L x 12 3/8 W x 1 1/4 H (in)	
CS-5001	




XL SIZE 8" PERFORATED BOTTOM	
Outside with Lid 30 3/8 L x 12 3/8 W x 7 2/3 H (in)	Usable Inside 28 1/2 L x 11 W x 5 1/2 H (in)
CS-5002	



XL SIZE CONTAINER BASKET	
Outside with Lid 28 1/8 L x 11 1/8 W x 6 1/4 H (in)	Usable Inside 28 L x 11 W x 5 3/8 H (in)
CS-5150	

# Special Features



HALF SIZE

THREE QUARTER

FULL

XL


ACCESSORIES

RESOURCES

**LONGEST CONTAINER FOR ROBOTIC INSTRUMENTS**


Fits robotic instruments without bending instrument tips

Container designed to hold multiple tray configurations: IN-8915-TF or IN-8916-TF



## ADDITIONAL FEATURES

- Combined weight does not exceed AMMI standard of 25 lbs
- Only available with perforated bottom



# Locks, Filters + Data Cards

HALF

THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES



## CONTAINER LOCKS WITH PROCESS INDICATOR

Quantity 1000
CS-0100



## CONTAINER FILTERS WITH PROCESS INDICATOR

Quantity 500
CS-0250



## DATA CARDS

Quantity 250
CS-0310



# Plaque Options



*Our in-house laser engraving capabilities allow us to personalize plaques with custom identification and/or specifications, such as interior instruments, exterior procedure types, specialty or facility information. Engrave up to 13 characters on a single line of text!*



HALF SIZE

THREE QUARTER

FULL

XL




ACCESSORIES

RESOURCES

## BLANK CONTAINER PLAQUES

### Color Options

Indicate plaque color using the corresponding numbers below in place of XX.

02	RED	
03	GREEN	
04	BLUE	
05	BLACK	
06	YELLOW	
07	VIOLET	
08	ORANGE	
09	BROWN	

CS-5500-XX

## ENGRAVED CONTAINER PLAQUES

### Color Options

Indicate plaque color using the corresponding numbers below in place of XX.

02	RED	
03	GREEN	
04	BLUE	
05	BLACK	
06	YELLOW	
07	VIOLET	
08	ORANGE	
09	BROWN	


CS-5500-XX-EN



# Metric Measurement Table

- HALF
- THREE QUARTER
- FULL
- XL
- ACCESSORIES
- RESOURCES

Container Size	Container Part Number	Outside with Lid (mm)	Usable Inside (mm)
HALF	CS-1001	299 L x 278 W x 32 H	N/A
	CS-1002	299 L x 278 W x 105 H	260 L x 260 W x 83 H
	CS-1003	299 L x 278 W x 143 H	260 L x 260 W x 114 H
	CS-1012	299 L x 278 W x 105 H	260 L x 260 W x 67 H
	CS-1013	299 L x 278 W x 143 H	260 L x 260 W x 102 H
	CS-1120	216 L x 248 W x 51 H	213 L x 244 W x 51 H
THREE QUARTER	CS-2001	464 L x 278 W x 32 H	N/A
	CS-2002	464 L x 278 W x 105 H	425 L x 260 W x 83 H
	CS-2003	464 L x 278 W x 143 H	425 L x 260 W x 114 H
	CS-2012	464 L x 278 W x 105 H	425 L x 260 W x 70 H
	CS-2013	464 L x 278 W x 143 H	425 L x 260 W x 102 H
	CS-2120	391 L x 248 W x 51 H	387 L x 244 W x 51 H
FULL	CS-3001	587 L x 278 W x 32 H	N/A
	CS-3002	587 L x 278 W x 105 H	549 L x 260 W x 83 H
	CS-3003	587 L x 278 W x 143 H	549 L x 260 W x 114 H
	CS-3004	587 L x 278 W x 208 H	549 L x 260 W x 181 H
	CS-3005	587 L x 278 W x 263 H	549 L x 260 W x 241 H
	CS-3012	587 L x 278 W x 105 H	549 L x 260 W x 70 H
	CS-3013	587 L x 278 W x 143 H	549 L x 260 W x 102 H
	CS-3014	587 L x 278 W x 208 H	549 L x 260 W x 168 H
	CS-3015	587 L x 278 W x 263 H	549 L x 260 W x 229 H
	CS-3120	502 L x 248 W x 52 H	498 L x 244 W x 51 H
	CS-3160	502 L x 248 W x 152 H	498 L x 244 W x 152 H
	XL	CS-5001	771 L x 315 W x 32 H
CS-5002		771 L x 315 W x 195 H	724 L x 279 W x 140 H
CS-5150		714 L x 283 W x 140 H	711 L x 279 W x 137 H



# Tray Cross Reference

Tray Part Number	Container Part Number	
	<i>SOLID</i>	<i>PERFORATED</i>
IN-1306-TF	CS-1002	CS-1012
IN-1310-TF	CS-2002	CS-2012
IN-1315-TF	CS-2002	CS-2012
IN-2006-TF	CS-1003	CS-1013
IN-2840-TF	CS-3004	CS-3014
IN-2841-TF	CS-3004	CS-3014
IN-2842-TF	CS-3004	CS-3014
IN-2847-TF	CS-3005	CS-3015
IN-2890-TF	CS-3004	CS-3014
IN-2950-TF	CS-1003	CS-1013
IN-5010-TF	CS-1003	CS-1013
IN-5011-TF	CS-2002	CS-2012
IN-6021-TF	CS-3003	CS-3013
IN-6022-TF	CS-3003	CS-3013
IN-6023-TF	CS-3003	CS-3013
IN-6024-TF	CS-3003	CS-3013
IN-6031-TF	CS-2003	CS-2013
IN-6500-TF	CS-3003	CS-3013
IN-7008-TF	CS-1002	CS-1012
IN-7010-TF	CS-1002	CS-1012
IN-7011-TF	CS-2003	CS-2013
IN-7014-TF	CS-2003	CS-2013
IN-7020-TF	CS-1002	CS-1012
IN-7030-TF	CS-1003	CS-1013
IN-7120-TF	CS-2002	CS-2012
IN-7123-TF	CS-2003	CS-2013
IN-7130-TF	CS-2003	CS-2013
IN-7140-TF	CS-2003	N/A
IN-7213-TF	CS-3002	CS-3012
IN-7220-TF	CS-3002	CS-3012
IN-7223-TF	CS-3003	CS-3013
IN-7230-TF	CS-3003	CS-3013
IN-7234-TF	N/A	CS-5002
IN-7240-TF	CS-3003	CS-3014
IN-7250-TF	CS-3004	CS-3014
IN-7260-TF	CS-3004	CS-3014
IN-7323-TF	N/A	CS-5002
IN-7344-TF	N/A	CS-5002
IN-7423-TF	N/A	CS-5002
IN-7430-TF	N/A	CS-5002
IN-7533-TF	N/A	CS-5002
IN-7582-TF	N/A	CS-5002

Tray Part Number	Container Part Number	
	<i>SOLID</i>	<i>PERFORATED</i>
IN-7583-TF	N/A	CS-5002
IN-7633-TF	N/A	CS-5002
IN-7644-TF	N/A	CS-5002
IN-8610-TF	CS-2002	CS-2012
IN-8613-TF	CS-2003	CS-2013
IN-8615-TF	CS-1002	CS-1012
IN-8616-TF	CS-1003	CS-1013
IN-8620-TF	CS-3003	CS-3013
IN-8621-TF	CS-3003	CS-3013
IN-8622-TF	CS-2003	CS-2013
IN-8623-TF	CS-2003	CS-2013
IN-8630-TF	CS-3003	CS-3013
IN-8632-TF	CS-3003	CS-3013
IN-8640-TF	CS-3003	CS-3013
IN-8642-TF	CS-3003	CS-3013
IN-8645-TF	CS-3003	CS-3013
IN-8721-TF	CS-2003	CS-2013
IN-8722-TF	CS-1002	CS-1012
IN-8750-TF	CS-3003	CS-3014
IN-8751-TF	CS-3003	CS-3013
IN-8761-TF	CS-3003	CS-3014
IN-8884-TF	CS-2002	CS-2012
IN-8889-TF	N/A	CS-5002
IN-8894-TF	N/A	CS-5002
IN-8897-TF	CS-1003	CS-1013
IN-8898-01-TF	N/A	CS-5002
IN-8899-TF	CS-2003	CS-2013
IN-8913-TF	CS-3003	CS-3014
IN-8914-TF	N/A	CS-5002
IN-8915-TF	N/A	CS-5002
IN-8916-TF	N/A	CS-5002
IN-8932-TF	N/A	CS-5002
IN-8936-TF	N/A	CS-5002
IN-8942-TF	N/A	CS-5002
IN-8943-TF	CS-1003	CS-1013
IN-8944-TF	N/A	CS-5002
IN-8983-TF	N/A	CS-5002
IN-8984-01-TF	N/A	CS-5002
IN-8986-TF	N/A	CS-5002
IN-8987-TF	N/A	CS-5002
IN-8988-TF	N/A	CS-5002
IN-8989-TF	N/A	CS-5002

HALF SIZE

THREE QUARTER

FULL

XL

ACCESSORIES

RESOURCES





















**InstruSafe®**,  
An Innovia Medical Brand

**FOR US + INTERNATIONAL ORDERS:**  
**Summit Medical LLC**  
**An Innovia Medical Company**  
815 Vikings Parkway  
Suite 100  
St. Paul, MN 55121  
USA

1.888.229.2875  
[www.instrusafe.com](http://www.instrusafe.com)

**FOR UK ORDERS:**  
**DTR Medical Ltd**  
**An Innovia Medical Company**  
17 Clarion Court, Enterprise Park  
Swansea SA6 8RF  
UK

+44 (0) 1792 797910  
[www.dtrmedical.com](http://www.dtrmedical.com)

---

InstruSafe Instrument Protection Trays are proudly made in St. Paul, Minnesota, USA.

InstruSafe® Sterilization Container Systems are proudly made in St. Paul, Minnesota, USA.

Summit Medical LLC is an FDA registered facility, ISO 13485:2016 certified and CMDCAS certified.

InstruSafe Trays are validated for multiple sterilization cycles to give our customers more options.

Sterilization Containers are validated for multiple sterilization cycles to give our customers more options.

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

da Vinci, SP, S, Xi, Single-Site and EndoWrist are registered trademarks of Intuitive Surgical, Inc. InstruSafe, Instru-Zyme, Hold-Its™, Hold-Downs™, Rounders Point Guards, and Guardian are trademarks and registered trademarks of Summit Medical.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.

All brands are trademarks or registered trademarks of their respective owners.

InstruSafe® is a registered trademarks of Summit Medical LLC.

©2023 Innovia Medical. All rights reserved. MM-IS-5006 | RevA | APRIL2024