

# FDA 510(K) Clearance

InstruSafe® Instrument Protection Trays are FDA cleared for each sterilization cycle listed below for use with wrap and/or containers as a sterile barrier.

Reference [instrusafe.com](http://instrusafe.com) for the most up to date information.






**INSTRUSAFE**

## STERILIZATION CYCLES + STERILE BARRIERS

STERILIZATION CYCLE	WRAP	AESCULAP® CONTAINER	GENESIS® CONTAINER
Steam 4 minute	✓	✓	✓
Ethylene Oxide	✓	✗	✓
STERRAD® 100S Standard	✓	✓	✗
STERRAD® 100NX® Standard	✓	✓	✗
STERRAD® 100NX® Express	✓	✗	✗
STERRAD® 100NX® Flex	✓	✓	✗
STERRAD® NX® Standard	✓	✓	✗
STERIS® AMSCO® V-PRO® 1, V-PRO® 1 Plus, V-PRO® maX + V-PRO® maX 2	✓	✓	✗

## STERILE BARRIERS + FOOT OPTIONS

Sterile barrier and storage needs determine the foot type built into a tray.

STERILE BARRIER	FOOT OPTION
<p><b>WRAP</b></p> 	<p>Rounders Corner Guards <b>-R</b></p> 
<p><b>CONTAINER</b></p> 	<p>Tall Feet <b>-TF</b></p> 

Summit Medical is an FDA registered facility, ISO 13485 : 2016 certified and CMDCAS certified. Genesis® is a registered trademark of CareFusion Corporation. Aesculap® is a registered trademark of B. Braun Company. STERRAD®, 100S®, 100NX® are registered trademarks of Advanced Sterilization Products (ASP) Division of Ethicon, Inc. STERIS®, AMSCO®, V-PRO® 1, V-PRO® 1 Plus, V-PRO® maX and V-PRO® maX 2 are registered trademarks of STERIS Corporation. InstruSafe® is a registered trademark of Summit Medical. ©2020 Summit Medical. All rights reserved.