


InstruSafe® Trays with Wrap + Containers



FDA 510(K) Clearance Quick Reference Chart

InstruSafe Instrument Protection Trays have been cleared for each sterilization cycle listed below when using wrap as a sterile barrier and/or as indicated, when using with rigid containers.

STERILIZATION CYCLE	WRAP	AESCALAP® RIGID CONTAINER	GENESIS® RIGID CONTAINER
Steam 4 minute	✓	✓	✓
Ethylene Oxide	✓	✗	✓
STERRAD® 100S Standard	✓	✓	✗
STERRAD® 100NX® Standard	✓	✓	✗
STERRAD® 100NX® Express	✓	✗	✗
STERRAD® 100NX® Flex	✓	✓	✗
STERRAD® NX® Standard	✓	✓	✗
STERIS® AMSCO® V-PRO® VPRO® 1 + V-PRO® maX	✓	✓	✗

All information is subject to change, please reference www.instrusafe.com for the most up to date information.

Genesis® is a registered trademark of CareFusion Corporation. Aesculap® is a registered trademark of B. Braun Company. STERRAD®, 100S®, 100NX® are registered trademarks of Advanced Sterilization Products (ASP) Division of Ethicon, Inc. STERIS®, V-PRO®, V-PRO® 1 and V-PRO® maX are registered trademarks of STERIS Corporation.

Summit Medical is an FDA registered facility, ISO 13485 : 2016 certified and CMDCAS certified.
 ©2019 Summit Medical. InstruSafe® is a registered trademark of Summit Medical. All rights reserved.

QUAL-4088 | RevF | AUG2019